INFORMATIONS IMPORTANTES

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COMPLICATIONS COURANTES

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Références

1. Anaplastic Large Cell Lymphoma (ALCL) In Women with Breast Implants: Preliminary FDA Findings and Analyses. Center for Devices and Radiological Health, U.S. Food and Drug Administration; January 2011 [cited 2017 26 February].
2. Anaplastic large cell lymphoma (ALCL). US Food and Drug Administration; 20 January 2016 [cited 2017 26 February].
3. Lymphome Anaplasique à Grandes Cellules associé aux implants mammaires (LAGC-AIM): Point sur les investigations en cours. Point d'Information.: L'Agence nationale de sécurité du médicament et des produits de santé (ANSM). 06 July 2016 [cited 2016 06 July].
4. Breast Implants: Expert Advisory Panel Advice on Association with Anaplastic Large Cell Lymphoma.: Australian Government Department of Health Therapeutic Goods Administration.; 20 December 2016 [cited 2017 26 February].
5. Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL): US Food and Drug Administration; 22 March 2017 [cited 2017 22 March].
6. Clemens MW, Medeiros LJ, Butler CE, Hunt KK, Fanale MA, Horwitz S, et al. Complete Surgical Excision Is Essential for the Management of Patients With Breast Implant-Associated Anaplastic Large-Cell Lymphoma. *J Clin Oncol*. 2016;34(2):160-8. Epub 2015/12/03. doi:10.1200/jco.2015.63.3412. PubMed PMID: 26628470; PubMed Central PMCID: PMCPCMC4872006 online at <http://www.jco.org>. Author contributions are found at the end of this article.
7. Clemens MW, Horwitz SM. NCCN Consensus Guidelines for the Diagnosis and Management of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Aesthetic surgery journal*. 2017. Epub 2017/02/12. doi: 10.1093/asj/sjw259. PubMed PMID: 28184418.
8. ASPS/ASAPS Update Breast Implant-Associated Anaplastic Large Cell Lymphoma(BIA-ALCL) Quick Facts and FAQs: The American Society For Aesthetic Plastic Surgery & American Society of Plastic Surgeons. Available from: <http://www.surgery.org/downloads/blasts/BIA-ALCL/>.
9. Breast implants - report cases of anaplastic large cell lymphoma (ALCL): Medicines and Healthcare products Regulatory Agency; 10 July 2014 [cited 2017 26 February].
10. Information on BIA-ALCL: American Society of Plastic Surgeons; 2017 [cited 2017 26 February].
11. Joint ASPS & ASAPS Statement on Breast Implant-Associated ALCL 2016 [cited 2017 26 February]. Available from: <http://www.surgery.org/downloads/private/joint-asps-asaps-statement-on-breast-implant-associated-alcl.pdf>.
12. Frequently Asked Questions (FAQ): A Guide to Breast Implant-Associated Anaplastic Large Cell Lymphoma: American Society of Plastic Surgeons and The Plastic Surgery Foundation; 2016 [cited 2017 22 March]. Available from: <http://www.thepsf.org/Documents/Clinical/PROFILE/profile-faq.pdf>.
13. International Society of Aesthetic Plastic Surgery (ISAPS) 2017. Available from: <http://www.isaps.org>.
14. MemoryShape Post-Approval Cohort Study (formerly Contour Profile Gel Core Study) Final Clinical Study Report. Mentor Worldwide, LLC; 02 June 2015.
15. MemoryGel Core Gel Clinical Study Final Report. Mentor Worldwide, LLC; April 2013.
16. Mentor MemoryShape Post-Approval Continued Access Study (formerly Contour Profile Gel Continued Access Study), Final Report. October 2014.
17. Mentor MemoryGel Breast Implant Large Post Approval Study Re-Op Phase Annual Report. 17 June 2016.
18. Adjunct Study Final Report for Mentor's MemoryGel Silicone Gel-filled Breast Implants. 02 November 2012.
19. Mentor MemoryShape CPG Styles Study: A Study of the Safety of the Contour Profile Gel Breast Implants in Subjects who are Undergoing Primary Breast Augmentation, Primary Breast Reconstruction, or Revision, Final Clinical Study Report. October 2015.
20. Clemens, MW. Breast Implant Associated-ALCL: US Experience on an Emerging Malignancy. International Expert Meeting on BIA-ALCL, Amsterdam, Netherlands. November 2018.
21. de Boer, M., et al., Breast implants and the risk of anaplastic large-cell lymphoma in the breast. *JAMA Oncology*, 2018. 4(3): p. 335-341
22. Brody, G.S., et al., Anaplastic Large Cell Lymphoma Occurring in Women with Breast Implants: Analysis of 173 Cases. *Plastic and Reconstructive Surgery*, 2015. 135(3): p. 695-705.
23. Gidengil, C.A., et al., Breast Implant-Associated Anaplastic Large Cell Lymphoma: A Systematic Review. *Plastic and Reconstructive Surgery*, 2015. 135(3): p. 713-720.
24. Loch-Wilkinson, A., et al., Breast implant associated Anaplastic Large Cell Lymphoma in Australia and New Zealand - high surface area textured implants are associated with increased risk. *Plastic and Reconstructive Surgery*, 2017. 140(4): p. 645-654.
25. Doren, E.L., et al., U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plastic and Reconstructive Surgery*, 2017. 139(5): p. 1042-1050.
26. Srinivasa, D.R., et al., Global Adverse Event Reports of Breast Implant-Associated ALCL: An International Review of 40 Government Authority Databases. *Plast Reconstr Surg*, 2017. 139(5): p. 1029-1039.
27. Johnson, L., et al., Breast implant associated anaplastic large cell lymphoma: The UK experience. Recommendations on its management and implications for informed consent. *Eur J Surg Oncol*, 2017. 43(8): p. 1393-1401.
28. Deva, A.K. "BIA-ALCL: Translating Science Into Practice." The Aesthetic Meeting of ASPAS, April 29, 2018, Javits Center, New York, NY. Lecture in Panel: Hot Topics in Breast Surgery-ALCL, Texture, Biofilms
29. Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL). US Food & Drug Administration. March 2018 [accessed 12DEC2018]. <https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/BreastImplants/ucm239995.htm>
30. HEALTH CANADA website information on ALCL <http://healthycanadians.gc.ca/recall-alert-rappel-avis/hc-sc/2017/65188a-eng.php>

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© Mentor Worldwide LLC 2019 CAN 109421-190424
Une unité de Produits médicaux Johnson & Johnson Inc.
Division de Johnson & Johnson Inc.