FMS[®]
Fluid Management System

5-Way Foot Pedal
**for FMS VUE[™] and NEXTRA[™] Fluid Management
and Tissue Debridement Systems**

Ref. 283573

P/N: 109724-KIM
Rev: C Revised 02/2017
© DePuy Synthes 2010-2017. All rights reserved.

FMS® 5-Way Foot Pedal

DESCRIPTION

The FMS 5-way Foot Pedal controls shaver and suction functions when used with the FMS VUE™ or NEXTRA™ Fluid Management Systems.

INTENDED USE

The FMS 5-Way Foot Pedal is intended for use with the FMS VUE or NEXTRA Fluid Management and Tissue Debridement Systems in a surgical setting by personnel trained in arthroscopy.

INDICATION FOR USE

The 5-Way Foot Pedal is an accessory to the FMS VUE and NEXTRA Fluid Management and Tissue Debridement Systems. The FMS VUE and NEXTRA Fluid Management and Tissue Debridement Systems which provide controlled fluid distention and suction, controlled cutting, burring, shaving and abrading of bone and tissue for arthroscopic surgery of the shoulder, knee, ankle, elbow, wrist and hip joints.

CONTRAINDICATION


Do not use the 5-Way Foot Pedal with systems other than the FMS VUE or NEXTRA Fluid Management Systems.

WARNINGS

- This equipment may be operated only by arthroscopic surgeons. The surgeon using this device must be trained in arthroscopic surgical procedures, be aware of the risks associated with those procedures, and have current knowledge of technological advances in surgical products and techniques.
- Prior to use, check the system components for damage. Check the cable integrity carefully. If there are signs of damage, do not use.
- Failure to follow all applicable instructions may result in serious surgical consequences to the patient.
- Factory settings should be used only as guidelines. The surgeon is responsible for settings pertaining to his surgical procedure.
- Do not wrap Foot pedal power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire and/or injury to patient or surgical personnel.
- This device has been tested and found in compliance with limits of medical devices according IEC 60601-1-2 standard. These limits are designed to provide reasonable protection in a typical medical installation, against harmful interference when devices are located in close proximity. If the pump is causing harmful interference with other devices, turn the main power off, then on again to determine if it is causing the interference, relocate the system or separate it from other devices. If you cannot resolve the problem contact DePuy Customer Service at 1-800-382-4682. In the European Union, contact your local affiliate.
- Disconnect the device from the pump when cleaning or inspecting.

PRECAUTIONS

- Electrical safety testing should be performed by a biomedical engineer or other qualified person.
- Read the instructions, cautions, and warnings provided with all FMS VUE System accessories before use. Your sales representative can advise which accessories are used with the FMS VUE and NEXTRA systems.

- 1 Reverse Pedal
- 2 Flow+™ Button
- 3 Oscillate Pedal
- 4 BLOODSTOP™ Button
- 5 Forward Pedal

INSTRUCTIONS FOR USE

Refer to the appropriate manual for instructions for use with the FMS VUE and NEXTRA Fluid Management and Tissue Debridement Systems.

FMS NEXTRA Instruction Manual - 109652

FMS VUE Fluid Management and Tissue Debridement System Operator's Manual - 110665

CLEANING AND DECONTAMINATION

After each operation, clean the foot pedal and the power cord with a Glutaraldehyde disinfectant such as CIDEX® OPA solution prepared and used as per manufacturer's instructions.

1. Prepare the disinfectant solution per manufacturer instructions.
2. Wipe the surface of the Foot Pedal with a cloth, swab or gauze pad soaked in disinfectant solution.
3. Let the disinfectant stay on all device areas for 19-20 minutes.
4. Rinse thoroughly in running water.
5. Allow the Foot Pedal to drain and dry thoroughly.

Caution:

DO NOT AUTOCLAVE

Do not twist the cable when storing. Do not unplug instrument by pulling on cable.

MAINTENANCE

Inspect all equipment and cables periodically for wear. If damage is noted, replace or return to Repair Service Center.

Any repairs must be performed by a qualified Biomedical Engineer.

STORAGE

Room temperature: -10°C (14°F) to +50°C (122°F)

Relative humidity: 10% to 90% condensation included.

Atmospheric pressure: between 500 hPa & 1060 hPa.

OPERATING CONDITIONS

Temperatures: +10°C (50°F) to +40°C (104°F)*

Relative humidity: 30% to 75%

Atmospheric pressure: between 700 hPa & 1060 hPa

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DePuy Mitek, Inc.
**DePuy International Ltd.**