

Fixed Bearing




 **Attune**[®]
Knee System

Revision Solutions

*stability***in***motion*[™]


BRINGING PATENTED TECHNOLOGIES TO A SEAMLESS SYSTEM, FROM PRIMARY THROUGH TO REVISION

The ATTUNE® Revision Fixed Bearing Knee System is a comprehensive system that is designed to enable you to effectively manage a broad range of complex primary and revision knee procedures.



ATTUNE
Primary Knee System

Breakthrough discoveries and proprietary technologies, clinically designed to address patient needs for stability and motion¹



The Revision Femoral Component is compatible with the ATTUNE Revision Fixed Bearing Tibial Insert and the ATTUNE Posterior Stabilized (PS) Fixed Bearing Insert.

ATTUNE
Revision Fixed Bearing Knee System

Comprehensive range of solutions that are designed to address the main challenges presented in complex primary and revision surgery: fixation², instability², and OR efficiency^{1,2}

ADDRESSING FIXATION AND PATIENT FIT





Bone Defect Compensation

Femoral Metaphyseal Sleeves feature a proprietary stepped design to compensate for substantial bone defects, compressively load the bone and provide a solid foundation for implant fixation.

360° Offset Capability

Offset options of 2, 4, and 6 mm with orientation available from 0° to 360° are designed to increase patient fit in both femoral and tibial anatomical variations.

Rotational Freedom

The ATTUNE Revision Cam/Spine interface is designed to provide +/- 1.25° varus/valgus constraint while allowing +/- 4° of internal/external rotational freedom during the full range of motion.

Improved Stem Geometry

The Press-Fit Stems provide rotational stability and create a balance between stiffness and flexibility designed to reduce stress associated with end-of-stem pain.³

Closing The Flexion Gap


The ATTUNE Revision Femoral Boss is designed to best fit a patient's natural anatomy through optimal placement in terms of flexion gap balance, bone coverage and valgus angle. The ATTUNE Revision Femoral Component has 1 mm additional posterior thickness, the same as the ATTUNE PS Femoral Component, to help close the flexion gap.

REDUCING INSTABILITY THROUGH IMPROVED KINEMATICS

Femoral-to-Insert Size Matching

The LOGICLOCK™ Tibial Base enables the Femoral Component and Tibial Insert to match size-to-size every time, allowing for optimization of the tibiofemoral contact mechanics and stability throughout the range of motion. Additionally, it is designed to accommodate 2-up and 2-down sizing between the Tibial Insert and Tibial Base.^{4,5,6}





Enhanced Tibiofemoral Conformity

ATTUNE GRADIUS™ Curve technology is designed to allow for enhanced tibiofemoral conformity where the gradually reducing femoral radius provides a smooth transition through the gait cycle.^{4,5,7}

Femoral Offset

360° offset capability with 2, 4, and 6 mm options on the Femoral Component developed to enable balancing the flexion gap while providing fixation and improved patient fit.

Tibio-femoral Sagittal Conformity (Conformity Ratio - $R_{\text{femoral}}/R_{\text{insert}}$)⁷

	0°	30°	60°	90°	120°
ATTUNE Revision Femoral Component with ATTUNE Revision FB Tibial Insert	0.88	0.73	0.59	0.40	0.20
ATTUNE Revision Femoral Component with ATTUNE Primary PS FB Tibial Insert	0.88	0.73	0.59	0.40	0.20
ATTUNE Primary PS Femoral Component with ATTUNE Primary PS FB Insert	0.88	0.72	0.59	0.40	0.20
SIGMA® TC3 Femoral Component with TC3 FB Insert	0.75	0.75	0.47	0.47	0.47
SIGMA TC3 Femoral Component with PS FB Insert	0.75	0.75	0.47	0.47	0.47

OR EFFICIENCY WITH STREAMLINED WORKFLOWS



Cut Through Trials

Cut Through Trials are designed to allow for setting rotation and balancing flexion and extension gaps in real-time. Cut Through Trials also have the ability to make Femoral Augment and Box resections and assemble the Augment and Box Trials without being removed.



Real-time Gap Assessment with Conventional Cut Guide

Conventional Cut Guides allow for assessing the anterior resection level, flexion gap, and M/L fit simultaneously by fine-tuning the femoral offset position in 360°. Additionally, Conventional Cut Guides replicate the distal thickness of the femoral implant, enabling assessment of extension gap before making bone resections.



360° Tibia Offset Workflow

The ATTUNE Revision Offset Tibial preparation allows for the Tibial Base to be positioned 360° around the Reamer to efficiently address anatomical variation.



Consistent Sleeve Preparation

Ream, broach, cut. Consistent Sleeve preparation promotes technique efficiency. Each Femoral Sleeve distalizes the joint line by 4 mm, matching the size increments of the Femoral Distal Augments.



Smart Kitting

Instrumentation kitted for streamlined surgical workflows while enabling management of a broad range of complex primary and revision knee procedures.

CONTINUING TO BALANCE THE PATIENT NEED FOR FREEDOM OF MOVEMENT AND JOINT STABILITY

The ATTUNE Knee System is designed to provide patients with STABILITY IN MOTION™, offering a comprehensive range of implant solutions that address overall patient satisfaction. In a large multi-center study, the ATTUNE Primary Knee System demonstrated improved outcomes across a broad range of Patient Reported Outcome Measures (PROMs) compared to certain existing knee brands at one year minimum follow-up.¹ The same technologies that have helped deliver these outcomes in primary are also incorporated in the ATTUNE Revision System. These features are designed to help address the percentage of patients who are not satisfied after knee revision surgery, which ranges from 27% to 38%.⁸





ATTUNE **GRADIUS™** Curve

The patented ATTUNE **GRADIUS™** Curve is a gradually reducing femoral radius designed to provide a smooth transition from stability to rotational freedom through a patient's range of motion.^{4,5,7}

GLIDERIGHT™ Articulation

The **GLIDERIGHT™** Articulation encompasses a trochlear groove designed to accommodate patient variation and soft tissue interaction, and patella components designed to optimize patella tracking while maintaining bone coverage.⁹

SOFCAM™ Contact

The Revision Tibial Inserts and Revision Femoral Component have taken into account the advantages of **SOFCAM™** Contact while considering the requirements for constraint in the Revision Construct throughout the range of motion. The AP kinematics for the Revision Femoral Component on the Revision Fixed Bearing Insert are the same as the Revision Femoral Component on the Primary PS Insert.¹⁰ The controlled engagement of the Cam to the Insert Spine provides a smooth transition from condylar control to Cam and Spine control.

LOGICLOCK™ Tibial Base

The **LOGICLOCK™** Tibial Base has a patented central locking design that provides the architecture for the system to optimize kinematics, while reducing backside micromotion to the lowest reported levels in the industry.¹¹

References:

1. Hamilton, W.G., Brenkel, I., Gibbon, A., Kantor, S., Clatworthy, M., Dwyer, K., Himden, S., Lesko, J. (2017). Early outcomes with a new primary total knee arthroplasty (TKA) system vs. contemporary TKA: Interim results of two worldwide, multi-center prospective studies. *American Academy of Orthopaedic Surgeons*, Poster 106.
2. Schroer, W.C., Berend, K.R., Lombardi, A.V., Barnes, C.L., Bolognesi, M.P., Berend, M.E., Ritter, M.A., Nunley, R.M. (2013). Why are total knees failing today? Etiology of total knee revision in 2010 and 2011. *The Journal of Arthroplasty*, 28 Suppl. 1: 116-119.
3. Leszko F., Gohsh, U., Heldreth, M., Barrett, D. (2014). The Effects of Different Revision TKA Stem Design on Bone Stress Distribution: FEA Comparative Study. *7th World Congress of Biomechanics*, Poster F143.
4. Clary, C.W., Fitzpatrick, C.K., Maletsky, L.P., Rullkoetter, P.J. (2012a, February). Improving dynamic mid-stance stability: an experimental and finite element study. *ORS Annual Meeting*, Poster #1044.
5. Clary, C.W., Fitzpatrick, C.K., Maletsky, L.P., Rullkoetter, P.J. (2013). The influence of total knee arthroplasty geometry on mid-flexion stability: an experimental and finite element study. *Journal of Biomechanics*, 46:1351-1357.
6. Fitzpatrick, C.K., Clary, C.W., Rullkoetter, P.J. (2012a). Post-cam engagement during dynamic activity with current posterior-stabilized TKR. *18th Congress of the European Society of Biomechanics (ESB)*, 1700: 29.
7. Fitzpatrick, C.K., Clary, C.W., Rullkoetter, P.J. (2012b, February). The influence of design on TKR mechanics during activities of daily living. *ORS Annual Meeting*, Poster #2034.
8. Greidanus, N.V., Peterson, R.C., Masri, B.A., Garbuz, D.S. (2011). Quality of life outcomes in revision versus primary total knee arthroplasty. *The Journal of Arthroplasty*, 26(4): 615-620.
9. Clary, C.W., Wright, A.P., Komosa, M.C., Maletsky, L.P. (2012b). The effect of patella medialization on patellofemoral kinematics after total knee replacement. *18th Congress of the European Society of Biomechanics (ESB)*, 29: 1262.
10. Shalhoub, S., Fitzwater, F., Dickinson, M., Clary, C.W., Maletsky, L.P. (2016). Quantifying the change in tibiofemoral kinematics between primary and revision total knee arthroplasty inserts. *International Society for Technology in Arthroplasty (ISTA) Annual Meeting*.
11. Leisinger, S., Hazebrouck, S., Deffenbaugh, D., Heldreth, M. (2011). Advanced fixed bearing TKA locking mechanism minimizes backside micromotion. *International Society for Technology in Arthroplasty (ISTA) Annual Meeting*.
12. Tokarski, Anthony T., et al. "Medicare fails to compensate additional surgical time and effort associated with revision arthroplasty." *The Journal of Arthroplasty*, 30.4 (2015): 535-538.



DePuy Orthopaedics, Inc.

700 Orthopaedic Drive
Warsaw, IN 46582
USA
Tel: +1 (800) 366-8143
Fax: +1 (800) 669-2530

DePuy International, Ltd.

St Anthony's Road
Leeds LS11 8DT
England
Tel: +44 (0) 113 270 0461

DePuy (Ireland)

Loughbeg, Ringaskiddy
Co. Cork, Ireland
Tel: +35 (321) 491 4000
Fax: +35 (321) 491 4199

Johnson & Johnson Medical Pty Ltd

t/a DePuy Synthes
1-5 Khartoum Road
North Ryde, NSW 2113
Australia

Johnson & Johnson (NZ) Ltd

507 Mt Wellington Road
Mt Wellington, Auckland 1060
New Zealand