DePuy Synthes

MITEK SPORTS MEDICINE

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

OMNISPAN™ Meniscal Repair System

C € 2797

P/N: 109282

Rev: E Revised 12/2019

© DePuy Synthes 2010-2019. All rights reserved.

ENGLISH

OMNISPAN™ Meniscal Repair System

CONTENTS

Contents are STERILE unless packaging is damaged or opened. Resterilization by any method is not recommended.

DESCRIPTION

The OMNISPAN™ Meniscal Repair System consists of the OMNISPAN backstop implants and needle, a sterile, disposable Deployment Gun, Malleable Graft Retractor, and the Arthroscopic Pusher/Cutter. The OMNISPAN backstop implant is a combination of two molded Polyetheretherketone (PEEK) implants, combined with ORTHOCORD® violet braided composite suture, Size 2/0. The molded implants together with the suture provide compression across the tear in the meniscus. The Deployment Gun properly introduces the implant into the meniscus and the Pusher/Cutter facilitates the final suture position flush with the meniscal surface.

SYSTEM MATERIALS

Backstops: Polyetheretherketone (PEEK)

Suture: ORTHOCORD violet braided composite suture, Size 2/0. ORTHOCORD Suture is a synthetic, sterile, braided, composite suture composed of dyed absorbable polydioxanone (PDS®) and undyed non-absorbable polyethylene. The partially absorbable suture is coated with copolymer of 90% caprolactone and 10% glycolide. The PDS copolymer has been shown to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption.

Needles: Stainless Steel with Silicone Tubing

Meniscal Deployment Gun: Polycarbonate, ABS, Ultem, and Stainless Steel (sold separately – covered by this IFU)

Malleable Graft Retractor: Stainless steel sold separately with the Meniscal Deployment gun – covered by this IFU).

Arthroscopic Pusher/Cutter: Stainless steel and aluminum (sold separately).

Note: The Arthroscopic Pusher/Cutter is provided non-sterile. Cleaning and sterilization instructions are covered by its own IFU that is packaged with it.

INDICATIONS

The OMNISPAN System is intended for use in the arthroscopic fixation of soft tissue procedures such as meniscal repair.

CONTRAINDICATIONS

1. Surgical procedures other than those listed in the INDICATIONS section.
2. Pathological conditions in the tissue which would impair secure fixation by suture.
3. Physical conditions which would eliminate, or tend to eliminate, adequate tissue strength, or retard tissue healing, i.e., blood supply limitation, infection.
4. Conditions which tend to preempt patient's ability to heal or the healing period, such as senility, mental illness or alcoholism are contraindicated.
5. Meniscal tears not suitable for repair because of the degree of damage to the meniscus body.
6. Comminuted bone surface, which would mitigate against secure fixation for the OMNISPAN implants, is contraindicated.
7. The OMNISPAN Meniscal Repair System is not designed for and should never be used to attach artificial ligaments.

WARNINGS

1. Users should be familiar with the arthroscopic surgical procedures and techniques for repair of meniscal tissue before using the OMNISPAN Meniscal Repair System.
2. Injury could result if device is not positioned correctly.
3. The OMNISPAN Meniscal Repair System (including backstop implants, delivery needle, ORTHOCORD suture, Deployment Gun and Malleable Graft Retractor) is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion, dulled edges; deformation and material degradation, which may impact the strength of the device and compromise device performance. Reprocessing of single-use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety. Discard opened and unused implants, Deployment Gun/malleable Graft Retractor.

PRECAUTIONS

1. A surgeon should not begin clinical use of the OMNISPAN Meniscal Repair System without reviewing the instructions for use and practicing the procedure in a skills laboratory.
2. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

3. Discard used needles in sharps container.
4. Use caution when tensioning the suture. Over-tensioning may cause suture or implant breakage.

MRI STATEMENT

The OMNISPAN™ Meniscal Repair System is MR Safe.

INSTRUCTIONS FOR USE

Assemble Needle to Deployment Gun

1. Open the implant package using aseptic technique.
2. Remove the needle from the package.
3. Pull up the Needle Lock Lever on the Deployment Gun to expose the chuck.
4. Insert the Deployment Gun shaft into the open end of the needle package and into the connector end of the needle until the needle CLICKS onto the Deployment Gun. Needle should always be attached to the Deployment Gun with the open slot in the needle facing upward.
5. Push down on the Needle Lock Lever to advance the sleeve over the needle and secure it in place.

Insertion of Implant into the Meniscus

1. Using a calibrated probe, measure the width of the meniscal tissue to be repaired.
2. Use a Malleable Graft Retractor to protect the implants and suture from getting caught on soft tissue, including the fat pad, during insertion into the knee. Once inside the joint space, remove the Malleable Graft Retractor. Alternatively, an arthroscopic cannula may be utilized.
3. Penetrate the tissue to desired depth using needle laser markings as a depth indicator (markings at 10mm, 15mm (double line), and 20mm).
4. At desired depth, squeeze the deployment lever (dark gray) while maintaining depth positioning to deliver the first implant. NOTE: There may be a slight pushback on the deployment gun while the implant goes through the tissue.
5. Remove the needle from the tissue and retain visibility of the tip of the needle in the scope to ensure the second implant maintains the proper position.
6. Pull the loading trigger (red) to load the second implant to the firing position of the needle (should be seen at 20mm marking).
7. Penetrate the tissue to desired depth using needle laser markings as a depth indicator.
8. At desired depth, squeeze the deployment lever (dark gray) while maintaining depth positioning to deliver the second implant. NOTE: There may be a slight pushback on the deployment gun while the implant goes through the tissue.

9. Release the deployment lever and remove the device from the joint space.
10. Pull the free suture leg with continuous force and a smooth motion until both suture bridge legs lay flat and tight against the tissue surface. NOTE: Minimize any starts and stops for optimal tensioning results.
11. Continue to implant devices to stabilize the tear as determined by the tear size following the above steps.
12. Once implants are placed, use the Arthroscopic Pusher/Cutter to cut the free end of the suture. Thread the suture through the instrument, advance the Arthroscopic Pusher/Cutter until flush with the repaired meniscus, maintain moderate tension on the suture limb and press the metal tab in the center ring of the Arthroscopic Pusher/Cutter to cut the suture. When using an instrument other than the Arthroscopic Pusher/Cutter to cut the suture, leave a limb of suture approximately 2mm.

ADVERSE EFFECTS

Adverse effects of implanted polymer devices may include inflammatory response and foreign body reactions.

STERILIZATION

The OMNISPAN Implant Repair System and Deployment Gun/Malleable Graft Retractor are provided sterile. The OMNISPAN implant Repair System is sterilized via EO method. The OMNISPAN Deployment Gun/Malleable Graft Retractor is sterilized via Gamma method. Do not resterilize.

IMPLANT STORAGE

Store in a cool dry place, below 25°C (77°F), away from moisture and direct heat. Once opened, the device should be used in surgery or discarded. Never return to storage. Do not use after expiration date.

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DePuy International Ltd.