**DePuy Synthes**

MITEK SPORTS MEDICINE

---

COMPANIES OF *Johnson & Johnson*

**SPEEDTRAP™**  
**Graft Preparation**  
**System**

**C € 2797**

**P/N: 113684**

**Rev: D Revised 12/2019**

© DePuy Synthes 2016-2019. All rights reserved.

# ENGLISH

## SPEEDTRAP™ Graft Preparation System

### DESCRIPTION

The SPEEDTRAP™ suture construct provides a means of preparing soft tissue grafts used in orthopedic procedures. The implantable suture construct is made of a non-absorbable, flat spine suture composed of ultra-high molecular weight polyethylene (UHMWPE) and a co-braided winding suture composed of either green (D&C Green #6) poly (ethylene terephthalate) (PET) and UHMWPE, or un-dyed PET and UHMWPE. The winding suture is coated with a non-absorbable, polybutylate coating. The SPEEDTRAP winding sutures meet USP requirements for Size 2 suture except for diameter. The suture construct is supplied on a disposable, non-implantable, sterile delivery device made of polypropylene.

## MATERIALS

- **Spine Suture:** UHMWPE
- **Winding Suture:** UHMWPE & Green PET (D&C #6), UHMWPE & un-dyed PET
- **Winding Suture Coating:** poly {oxy-1,4 butanediyl oxy (1,6-dioxo-1,6 hexanediyl)} or Polybutilate
- **Delivery Device:** Polypropylene

## INDICATIONS

The SPEEDTRAP suture construct is indicated for use in soft tissue approximation in orthopedic procedures.

## STERILIZATION

Contents are sterile unless package has been opened or damaged. Discard opened packages and un-used devices and suture constructs.

## CONTRAINDICATIONS

- SPEEDTRAP Graft Preparation System is not indicated for primary fixation in distal biceps procedures.
- SPEEDTRAP Graft Preparation System is not indicated for primary fixation procedures where the sutures are not kept under tension.
- A single SPEEDTRAP Graft Preparation System should not be used to join two tendon ends together.

## **ACTIONS**

The SPEEDTRAP suture is a permanent, non-absorbable construct which may become encapsulated in the surrounding connective tissue.

## **WARNINGS**

- Do not trim graft ends less than 3mm from the distal whip stitch; the suture construct may slip off the graft end.
- Do not use the SPEEDTRAP Graft Preparation System on a graft larger than 12mm, or on combined grafts with a total diameter larger than 12mm.
- This product is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.

## **PRECAUTIONS**

- Deploy the SPEEDTRAP Graft Preparation delivery device external to the surgical site.
- Orient the SPEEDTRAP Graft Preparation delivery device with the suture tails distal to the graft end; otherwise, difficulties with passing the graft through tunnels may result.
- Do not manipulate sutures prior to use.

- Post-surgery, and under certain circumstances, immobilization by external support may be employed at the discretion of the surgeon.
- Do not use the SPEEDTRAP suture construct on individual graft diameters smaller than 2mm.
- Avoid damage to the device or sutures when handling. Avoid crushing or crimping damage caused by application of surgical instruments, such as forceps or graft preparation board clamps.
- As with other suture materials, adequate knot security requires the standard surgical technique of flat and square ties with additional throws if indicated by surgical circumstance and the experience of the operator.
- As with other suture tensioning techniques, avoid under-tensioning the SPEEDTRAP winding sutures; otherwise, laxity of the SPEEDTRAP suture/graft construct may result. Similarly, avoid over-tensioning as this may compromise the integrity of the SPEEDTRAP suture/graft construct.

## **ADVERSE REACTIONS**

Adverse effects associated with the use of this type of suture construct may include infection and minimal acute inflammatory tissue reaction.

## **MRI STATEMENT**

The SPEEDTRAP Graft Preparation System is MR Safe.

## **INSTRUCTIONS FOR USE**

- Orient the delivery device with winding suture tails distal to the graft end. Use the same orientation when preparing a graft end or a folded portion of a graft.
- Pass the graft through all suture loops until the graft end is flush with jaws of the delivery device.
- Clamp the delivery device around graft and hold the clamped position.
- Apply axial tension to the graft while pulling each winding suture limb until all suture loops are tight around the graft.
- Unclamp the delivery device, remove from graft and discard delivery device.

## **HOW SUPPLIED**

The SPEEDTRAP suture construct is supplied sterile, braided and colored either Green/White-striped or White strands, Size 2 (metric 5) in single or multi-packs in boxes of 4.

## **STORAGE**

Sterile product, once opened, use in surgery or discard. Never return to storage.

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