

GLOBAL™ APG GLOBAL™ STEPTECH®

ANCHOR PEG GLENOID

Design Rationale



GLOBAL™ APG

SHOULDER SYSTEM

The modular shoulder system that provides surgeons adaptability in the Operating Room without compromise.



GLOBAL UNITE™
Platform Shoulder System



GLOBAL ICON®
Stemless Shoulder System

GLOBAL™ APG GLOBAL™ STEPTECH® ANCHOR PEG GLENOID

The GLOBAL™ Anchor Peg Glenoid Implant System and the GLOBAL™ STEPTECH are part of the DePuy Synthes GLOBAL Implant Family. They are compatible with both GLOBAL UNITE™ and the stemless GLOBAL ICON®.



GLOBAL™ APG ANCHOR PEG GLENOID

The patented GLOBAL Anchor Peg Glenoid Implant System is addressing the need for long-term implant fixation and stability, based on a proven design:

100%

survivorship at an average of 34 months¹

69%

reduced risk of revision with GLOBAL APG at 11 years²

4%

GLOBAL APG cumulative percentage revision at 10 years versus 13.4% all other stemmed shoulders²






Less primaries revised due to loosening with GLOBAL APG compared to all other stemmed shoulders²



Addressing glenoid fixation

Proven design

- Fluted central peg facilitates bony integration around the flutes³
- Cemented pegs are able to lessen glenoid loosening and resist shear forces at the bone-cement-implant interface⁴
- Three cemented peripheral pegs are designed to provide immediate stability in the bone³
- 5 sizes ranging from 40-56 in 4 mm increments



Advanced biomechanics

The APG glenoid is designed to create a 6 mm mismatch between the glenoid and the humeral head component.⁵ As demonstrated in a cadaveric study, this promotes a more natural range of motion by emulating the anatomic biomechanics of a shoulder.⁶ Mismatched designs have shown to have less surface damage, less radiolucency and less impingement and edge abrasion compared to conforming designs.⁷

Advantages of all-polyethylene PREMIERON X-Linked Polyethylene glenoids

- Decreased rate of particle generation and wear-particle induced osteolysis⁴
- Potential for improved implant survival and clinical outcomes⁴
- Optimized material property for multidirectional wear, improved wear performance⁸

GLOBAL™ STEPTECH®

ANCHOR PEG GLENOID



Addressing glenoid erosion in TSA

The stepped, augmented glenoid is a modification of the all-polyethylene pegged GLOBAL™ Anchor Peg Glenoid featuring the cemented peripheral pegs and the fluted central peg, designed to provide long-term fixation through an interference fit and bony ingrowth. It is designed to compensate for posterior glenoid erosion or bone deficiency, removing the need for bone grafting and reducing the need for asymmetric reaming, hence streamlining the surgical procedure. Patients who have Type B1 or B2 glenoids may benefit most from this device.


Posterior glenoid bone loss is a common finding in advanced glenohumeral arthritis and failure to address this bone loss during TSA will likely lead to suboptimal results. Re-establishing ideal glenoid version during TSA can be challenging in presence of significant posterior glenoid wear and literature suggests retroversion greater than 10° should be corrected.⁹ Glenoid bone loss and increased retroversion are significantly associated with glenoid component loosening.

Addressing the clinical need

Posterior Glenoid Erosion – B2 Glenoid

Asymmetric wear can lead to increased reduction of the posterior surface of the glenoid fossa. Studies show asymmetric wear can result in increased retroversion in osteoarthritic patients.¹⁰

Augmented glenoid components offer a solution to this difficult pathology. All-polyethylene augmented glenoid components with a posterior step performed favorably compared to traditional methods (reaming or grafting).¹¹



Product Overview



Step Height	Version Correction (approximate)
+3 mm	5°
+5 mm	10°
+7 mm	15°


- Designed to restore glenoid version and to offset posterior glenoid bone loss¹¹ without need for medialization, reaming and implant under sizing or bone grafting
- The amount of correction is managed by the heights (+3, +5 or +7 mm) of the posterior step built into the component.

- Supports the need to place the component in an optimal version position
- Same peg fixation design and implant range (40-56 mm) as the APG
- These features work together to decrease the shear stresses at the bone/implant and bone/cement interface, thereby increasing implant stability as demonstrated by computational modeling.¹¹

Addressing anterior glenoid lift-off



Different augmented glenoid designs vary in their resistance to loosening in mechanical testing. Results from these tests show that Steptech is the only augmented glenoid design that does not have a higher degree of lift-off than the standard APG for the central peg after performance testing, which may make it more likely to maintain long-term fixation.

As seen in the image, the posterior step is 77° to perpendicular. As shown in mechanical testing, this design counteracts posterior loading.¹² Augmented component with a surface that articulates with the prepared bone surface to be perpendicular to the vector of joint loading provides better mechanical properties.¹²



TRUMATCH[®]

PERSONALIZED SOLUTIONS SHOULDER SYSTEM



The TRUMATCH Personalized Solutions Shoulder System, enables surgeons to analyze the patient's glenohumeral anatomy in multiple dimensions and perform pre-surgical planning in a 3D environment. The 3D Reconstruction of the glenoid can be used to help plan for which implant step height and amount of needed correction will be needed.



References:

1. Groh GI, Survival and radiographic analysis of a glenoid component with a cementless fluted central peg. J Shoulder Elbow Surgery (2010) 19: 1265-1268.
2. AOA Automated Industry Report 3608 Depuy Synthes Australia Global Total Stemmed Shoulder. Report Generated: 9 February 2021. Data Period: 16 April 2004 – 05 Februar 2021.
3. Wirth MA, et al. Total shoulder arthroplasty with an all polyethylene pegged bone-ingrowth glenoid component: a clinical and radiographic outcome study. J Bone Joint Surg Am. 2012 Feb 1;94(3):260-7.
4. Carpenter SR, Urits I, Murthi AM. Porous metals and alternate bearing surfaces in shoulder arthroplasty. Curr Rev Musculoskelet Med. 2016 Mar;9(1):59-66.
5. IFU-0902-00-457 Rev N
6. Karduna AR, Williams GR, Williams JL, Iannotti JP. Glenohumeral joint translations before and after total shoulder arthroplasty. A atudy in cadavera. Journal Bone Joint Surgery Am 1997;79:1166-74
7. Nho SJ, et al. Comparison of conforming and nonconforming retrieved glenoid components. J Shoulder Elbow Surg. 2008 Nov- Dec;17(6):914-20.
8. Wirth MA, et al. Cross-linked glenoid prosthesis: A wear comparison to conventional glenoid prosthesis with wear particulate analysis. Journal of shoulder and elbow surgery. 2009; 18, 130-137.
9. Favorito PJ, et al. Total shoulder arthroplasty for glenohumeral arthritis associated with posterior glenoid bone loss: results of an all-polyethylene, posteriorly augmented glenoid component. Journal of shoulder and elbow surgery. 2016; 25, 1681-1689.
10. Walch G, Badet R, Boulahia A, Khoury A. J Arthroplasty 1999; 14(6):756-60.
11. Sabesan VJ, et al. Biomechanical comparison of 2 augmented glenoid designs: an integrated kinematic finite element analysis. Journal of shoulder and elbow surgery. 2019: 1166-1174.
12. Iannotti JP, et al. Liftoff resistance of augmented glenoid components during cyclic fatigue loading in the posterior-superior direction. Jpurnal of shoulder and elbow surgery. 2013; 1-7.

This publication is not intended for distribution outside the EMEA region.

The third-party trademarks used herein are the trademarks of their respective owners.

For full product details including indications, contraindications, warnings and precautions, please consult the Instructions for Use.



Johnson & Johnson Medical Limited. Baird House, 4 Lower Gilmore Bank, Edinburgh, EH3 9QP, United Kingdom. Incorporated and registered in Scotland under company number SC132162.

DePuy Orthopaedics, Inc.
700 Orthopaedic Drive
Warsaw, IN 46582
USA
Tel: +1 (800) 366 8143
Fax: +1 (800) 669 2530

DePuy International Ltd
St Anthony's Road
Leeds LS11 8DT
England
Tel: +44 (0)113 270 0461

DePuy Ireland UC
Loughbeg
Ringaskiddy
Co. Cork
Ireland
Tel: +353 21 4914 000
Fax: +353 21 4914 199

www.depuySynthes.com