EMBOVAC™ Aspiration Catheter

Medos International SÀRL
Chemin-Blanc 38, 2400
Le Locle, Switzerland

CE
2797

LCN-208913-001/A
© 2019 Cerenovus.
All rights reserved.
♻️ Printed 03/19

ENGLISH

IMPORTANT INFORMATION

Please read before use

EMBOVAC™ Aspiration Catheter

STERILE

EO

Description

The EMBOVAC™ Aspiration Catheter is a single lumen, variable stiffness catheter. The catheter has a hydrophilic coating to reduce friction during use. The catheter includes a radiopaque marker on the distal end for angiographic visualization and a luer hub on the proximal end allowing attachments for flushing and aspiration. A hemostasis valve and two peelable introducers are included in the package.

Indications

The EMBOVAC Aspiration Catheter is indicated for general intravascular use in the neuro vasculature. The catheter can be used to facilitate introduction of diagnostic or therapeutic agents and is also intended for use in removal/ aspiration of emboli and thrombi from selected blood vessels in the neuro vasculature.

Contraindications

None known.

WARNINGS

- The EMBOVAC catheter should only be used by physicians trained in interventional endovascular procedures at medical facilities with the proper imaging equipment.
- Limited testing has been performed with solutions such as contrast media and saline. The use of these catheters for delivery of solutions other than the types that have been tested for compatibility is not recommended.
- Do not use with Ethiodol or Lipiodol contrast media, or other such contrast media which incorporates the components of these agents.
- The appropriate anti-coagulation and anti-platelet therapy should be administered in accordance with standard medical practice.
- This device is not intended use for use with power injectors.
- Excessive aspiration with the distal tip of the EMBOVAC catheter covered by the vessel wall may cause vessel injury. Carefully investigate location of the distal tip under fluoroscopy prior to aspiration.
- If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter damage or patient injury. Remove and replace catheter.
- Do not advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the device against resistance could dislodge a clot, perforate a vessel wall, or damage the device.
- Torqueing the catheter excessively while kinked may damage the device, resulting in separation of the catheter shaft. Withdraw the entire device (the device, microcatheter, and guidewire) if the device is severely kinked.
- Dispose of all used devices in accordance with hospital policy for biohazardous materials.
- To prevent thrombus formation and contrast media crystal formation, maintain a constant infusion of appropriate flush solution through catheter lumen.

Precautions

Inspect the sterile package carefully. Do not use if:

- o The package or seal appears damaged,
- o Contents appear damaged, or
- o The expiry date has passed.

- Carefully inspect all devices prior to use. Verify size, length, and condition are suitable for the specific procedure. Do not use a device that has been damaged in any way; replace with another EMBOVAC catheter. A damaged device may cause complications.
- Exercise care in handling the EMBOVAC catheter to reduce the chance of accidental damage.
- To control the proper introduction, movement, positioning and removal of the catheter within the vascular system, users should employ standard clinical angiographic and fluoroscopic practices and techniques throughout the interventional procedure.
- If removed from the patient, the hydrophilic coating on the EMBOVAC catheter should be hydrated with heparinized saline. Do not allow the coating to dry.

Adverse Events

Potential adverse events associated with the use of catheters or with the endovascular procedures include, but are not limited to:

- Access site complications
- Allergic reaction
- Aneurysm perforation
- Aneurysm rupture
- Arrhythmia
- Cardiac complications
- Cerebral edema
- Continued or re-occlusion of targeted arterial segment
- Death
- Embolism (air, foreign body, plaque, thrombus)
- Hematoma
- Hemorrhage
- Hypo/hypertension
- Infection
- Ischemia
- Neurological deficits
- Pseudoaneurysm
- Renal insufficiency
- Respiratory complications
- Seizure
- Stroke
- Transient Ischemic Attack
- Vasospasm
- Vessel dissection
- Vessel occlusion
- Vessel perforation
- Vessel rupture
- Vessel thrombosis

How Supplied

The EMBOVAC Aspiration Catheter is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Cerenovus will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner pouch is not opened or damaged, the product is sterile and nonpyrogenic.

Preparations for use

1. Set-up continuous flush through guide sheath, balloon guide, or guide catheter lumen.
2. Select an appropriately sized catheter based on intended procedure and anatomy.

3. Before removing the EMBOVAC Catheter from the packaging hoop dispenser, rotate the luer on the packaging hoop 90 degrees. Attach a syringe containing minimum 20 cc heparinized saline to the luer connector at the end of the hoop dispenser. Flush the dispenser to hydrate the hydrophilic coating of the catheter. Failure to do so can compromise the coating and lubricity of the catheter.
4. Gently remove the catheter and accessories from the hoop and inspect prior to use to verify that they are undamaged.
CAUTION: Do not use a catheter that has been damaged in any way. If damage is detected, replace with another EMBOVAC catheter that is not damaged.
5. Flush the catheter lumen with heparinized saline solution.
6. Attach compatible hemostasis valve with continuous flush to the catheter.

Directions for Use

1. Insert an appropriately sized guidewire or compatible microcatheter and guidewire into the EMBOVAC catheter and advance until the guidewire and the catheter are aligned at the distal end.
2. Using the introducer sheath provided in the package, insert the catheter and guidewire through a hemostatic valve of the guide sheath, balloon guide, or guide catheter.
3. Remove the peel away introducer sheath from the hemostasis valve and peel the introducer from the catheter once the distal shaft of the EMBOVAC catheter is placed inside the patient body.
4. Under fluoroscopic guidance, advance the catheter through the vasculature to the desired location.
WARNING: Do not advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy.

Aspiration through the EMBOVAC catheter

WARNING: Torqueing the EMBOVAC catheter excessively while kinked may damage the device, resulting in separation of the catheter shaft. Withdraw the entire device (the device, microcatheter, and guidewire) if the device is severely kinked.

1. Under fluoroscopy, confirm that the distal tip of the EMBOVAC catheter is engaged with the emboli or thrombi.
2. Remove the perfusion line from the proximal end of the EMBOVAC catheter and attach a stopcock with 60 cc locking syringe or aspiration system to the sidearm of the hemostasis valve.
3. Tighten the hemostasis valve.
4. With the stopcock closed, pull and lock the plunger of the syringe or turn on the aspiration system.
5. Open the stopcock to begin aspiration. Ensure that the syringe/aspiration system aspirates blood, emboli, or thrombi through the system. Close the stopcock if the syringe/aspiration system does not aspirate any blood, emboli, or thrombi, or if slow aspiration is observed. Carefully investigate the cause of the restriction and reposition the distal tip if necessary. If blood, emboli, or thrombi is stuck in the EMBOVAC catheter, remove the entire device (the device, microcatheter, and guidewire) and clear the inner lumen. If restriction remains and the distal tip of the catheter has been repositioned correctly, close the stopcock, re-attach the syringe/aspiration system, and resume aspiration. Maintain aspiration so that the emboli or thrombi remains fully engaged with the distal tip of the catheter. With the emboli or thrombi fully engaged, pull back the EMBOVAC catheter out of the patient body.
WARNING: Excessive aspiration with the distal tip of the EMBOVAC catheter covered by the vessel wall may cause vessel injury.
WARNING: If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter damage or patient injury. Remove and replace catheter.
6. After aspiration is completed, remove the device from the patient body. If re-access to the vasculatures with the same device is desired, flush and clean the inner lumen of the device by infusion. Inspect the device for any damage.
CAUTION: Do not use the device if any damage or irregularities are observed.

The physician has the discretion to modify described manipulations of the EMBOVAC catheter to accommodate the complexity and variation in procedures. Any technique modification must be consistent with previously described instructions, warnings, precautions and patient safety information.

Storage

Store in a cool, dark, dry place.

Warranty

Cerenovus warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

TM EMBOVAC is a trademark of Cerenovus or affiliated companies.

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