

Part of the DePuy Synthes Locking Compression Plate (LCP®) System

Large Fragment LCP® Instrument and Implant Set

Surgical Technique



Table of Contents

Introduction	Large Fragment LCP Instrument and Implant Set	2
	Locking Screws	4
	AO Principles	5
	Indications	6
	Fixation Principles—Conventional Plating	7
	Fixation Principles—Locked Plating	9
	Combining Conventional and Locked Plating Techniques	10
	Screw Selection	11

Surgical Technique	Preparation and Reduction	12
	Screw Insertion	14
	Screw Placement Verification	17
	Postoperative Treatment	18
	Implant Removal	18

Product Information	Implants	19
	Instruments	21
	Set Lists	23
	Discontinued Product	29

MR Information

The Large Fragment LCP Implant Set has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. The safety of the Large Fragment LCP Implant Set in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

 Image intensifier control

Large Fragment LCP Instrument and Implant Set

The aim of any surgical fracture treatment is to reconstruct the anatomy and restore its function. According to the AO, internal fixation is distinguished by anatomic reduction, stable fixation, preservation of blood supply and early, functional mobilization. Plate and screw osteosynthesis has been established as clinically beneficial for quite some time. Clinical results have improved by using internal fixation with angular stability (internal fixators) in complicated fractures and in osteopenic bone.

The DePuy Synthes Locking Compression Plate (LCP®) System is part of a Stainless Steel and Titanium Plate and Screw System that merges locking screw technology with conventional plating techniques. The LCP System has many similarities to conventional plate fixation methods, but with a few important improvements. Locking screws provide the ability to create a fixed-angle construct while utilizing familiar AO plating techniques. A fixed-angle construct provides improved fixation in osteopenic bone or multifragment fractures where traditional screw purchase is compromised. LCP Plate constructs do not rely on plate-to-bone compression to maintain stability, but function similarly to multiple small multiangled blade plates.

The following points distinguish treatment using Locking Compression Plate technology:

- Allows fracture treatment using conventional plating with conventional cortex or cancellous bone screws;
- Allows fracture treatment using locked plating with bicortical or unicortical locking screws;
- Permits the combination of conventional and locking screw techniques.

Note: The LCP System applies to many different plate types and is therefore suitable for a large number of fracture types. For that reason, this technique guide does not deal with any specific fracture type. For more information, please refer to *AO Principles of Fracture Management*,¹ to AO courses (www.ao-asif.ch), and to the appropriate plate-specific technique guide.

1. Rüedi TP, et al, ed., *AO Principles of Fracture Management*, New York: Thieme. 2000.

Locking compression plates

The LCP Plates have these LC-DCP® Plate features:

- Uniform hole spacing
- Load (compression) and neutral screw positions

Combination locking and compression Combi holes

The Combi holes allow placement of conventional cortex and cancellous bone screws on one side or threaded conical locking screws on the opposite side of each hole.




- Threaded hole section for locking screws
- DCU hole section for conventional screws
- Locking screw in threaded side of Combi hole
- Cortex screw in compression side of Combi hole

Notes:

- **Combi holes in straight plates are oriented with the conventional portion of each hole further from the middle of the plate. This facilitates utilization of the LCP Plates for dynamic compression using traditional AO techniques.**
- **The holes in the straight LCP Plates are larger at the two ends to allow the insertion of cancellous bone screws.**




224.591



4.0 mm and 5.0 mm locking screws, self-tapping, with StarDrive™ Recess

The locking screws mate with the threaded portion of the Combi holes to form a fixed-angle construct.



Locking Screws

Locking screw design

The screw design has been modified, from standard 4.5 mm cortex screw design, to enhance fixation and facilitate the surgical procedure.

New features include:

Conical screwhead

The conical head facilitates alignment of the locking screw in the threaded plate hole to provide a fixed-angle connection between the screw and the plate.

Large core diameter

The large core diameter improves bending and shear strength of the screw, and distributes the load over a larger area in the bone.

Thread profile

The shallow thread profile of the locking screws results from the larger core diameter, and is acceptable because locking screws do not rely solely on screw purchase in the bone to maintain stability.


Drive mechanism

The StarDrive Recess provides improved torque transmission to the screw, while retaining the screw without the use of a holding sleeve.

Uncortical screw fixation

Bicortical screw fixation has long been the traditional method of compressing a plate to the bone where friction between the plate and the bone maintains stability. Screw stability and load transfer are accomplished at two points along the screw: the near and far cortices.

Uncortical locking screws provide stability and load transfer only at the near cortex due to the threaded connection between the plate and the screw. Screw stability and load transfer are accomplished at two points along the screw: the screwhead and near cortex. Because the screw is locked to the plate, fixation does not rely solely on the pullout strength of the screw or on maintaining friction between the plate and the bone.



- a. Bicortical screws require two (2) cortices to achieve stability
- b. Uncortical screws utilize the locked screwhead and the near cortex to achieve stability

AO Principles


In 1958, the AO formulated four basic principles, which have become the guidelines for internal fixation.^{1,2}

Anatomic reduction

Fracture reduction and fixation to restore anatomical relationships.

Early, active mobilization

Early and safe mobilization and rehabilitation of the injured part and the patient as a whole.



Stable fixation

Fracture fixation providing absolute or relative stability, as required by the patient, the injury, and the personality of the fracture.

Preservation of blood supply

Preservation of the blood supply to soft tissues and bone by gentle reduction techniques and careful handling.

1. Rüedi TP, RE Buckley, CG Moran. *AO Principles of Fracture Management*. 2nd ed. Stuttgart, New York: Thieme; 2007.

2. Müller ME, Allgöwer M, Schneider R, Willenegger H. *Manual of Internal Fixation*. 3rd ed. Berlin, Heidelberg, New York: Springer-Verlag; 1991.

Indications

The Synthes Locking Compression Plates—Narrow and broad, are intended for fixation of various long bones, such as the humerus, femur, and tibia. They are also for use in fixation of periprosthetic fractures, osteopenic bone, and nonunions or malunions.

The Synthes Locking Compression Plates—T-Plates are intended to buttress metaphyseal fractures of the proximal humerus, medial tibial plateau, and distal tibia. They are also for use in fixation of osteopenic bone and fixation of nonunions and malunions.

The Synthes 4.5 mm LCP Proximal Tibia Plate is intended for treatment of osteopenic bone, tibial osteotomies, nonunions, malunions, and fractures of the proximal tibia including:

- Simple, comminuted fractures
- Lateral wedge fractures
- Depression medial wedge fractures
- Bicondylar combination of lateral wedge and depression fractures
- Periprosthetic fractures
- Proximal fractures with associated shaft fractures.

The Synthes 3.5 mm/4.5 mm LCP Medial Proximal Tibia Plates are intended to buttress metaphyseal fractures of the medial tibia plateau, split-type fractures of the medial tibia plateau, medial split fractures with associated depressions, and split or depression fractures of the medial tibia plateau. The plates may also be used for fixation of the proximal quarter (lateral and medial) of the tibia as well as segmental fractures of the proximal tibia. The 4.5 mm version may also be used for fixation of nonunions and malunions of the medial proximal tibia and tibia shaft, as well as opening and closing wedge tibial osteotomies.


Fixation Principles—Conventional Plating

The following pages show the biomechanical features of conventional plating techniques, locked or bridge plating techniques, and a combination of both.

Note: Please refer also to the *AO Principles of Fracture Management*,¹ to AO courses (www.ao-asif.ch), and to the corresponding special literature.


Construct stability

The tensile force (f_1) originating from tightening the screws compresses the plate onto the bone (f_2). The developing frictional force (f_3) between the plate and the bone leads to stable plate fixation. To ensure construct stability, the frictional force must be greater than the patient load (f_4).



Anatomic contouring of the plate

The aim of internal fixation is anatomic reduction, particularly in articular fractures. Therefore, the plate must be contoured exactly to match the shape of the bone.




Lag screw

Interfragmentary compression is accomplished with a lag screw. This is particularly important in intra-articular fractures which require a precise reduction of the joint surfaces. Lag screws can be angled in the plate hole, allowing placement of the screw perpendicular to the fracture line.


Primary loss of reduction

In conventional plating, even though the bone fragments are correctly reduced prior to plate application, fracture dislocation will result if the plate does not precisely fit the bone. In addition, if the lag screw is not placed perpendicular to the fracture line, shear forces will be introduced. These forces may also cause loss of reduction.



Secondary loss of reduction


Under axial load, postoperative secondary loss of reduction may occur by toggling of the screws in the plate. Since cortex screws do not lock to the plate, the screws cannot oppose the acting force and may loosen, or be pushed axially through the plate holes.



1. Rüedi TP, et al, ed., *AO Principles of Fracture Management*, New York: Thieme. 2000.

Blood supply to the bone

Construct stability depends upon compressing the plate to the bone. Therefore, the periosteum is compressed under the plate, reducing or even interrupting blood supply to the bone. The result is delayed bone healing due to temporary osteoporosis underneath the plate.



Osteoporosis

Due to a compromised cortical structure, screws cannot be tightened sufficiently to obtain the compression needed for a stable construct. This may cause loosening of the screws and loss of stability, and may cause loss of reduction.

Conventional plating achieves good results in:

1. Good quality bone; and
2. Fractures which are traditionally fixed with lag screws to achieve direct bone healing.

Special attention must be paid to:


1. Osteoporotic bone—during rehabilitation, the load should be kept to a minimum to prevent postoperative loss of reduction; and
2. Multifragment fractures—the anatomic reduction may be accomplished at the expense of extensive soft tissue trauma and denudation.

Fixation Principles—Locked Plating

Screws lock to the plate, forming a fixed-angle construct.

Maintenance of primary reduction


Once the locking screws engage the plate, no further tightening is possible. Therefore, the implant locks the bone segments in their relative positions regardless of degree of reduction.



Precontouring the plate minimizes the gap between the plate and the bone, but an exact fit is not necessary for implant stability. This feature is especially advantageous in minimally or less invasive plating techniques because these techniques do not allow exact contouring of the plate to the bone surface.


Stability under load

By locking the screws to the plate, the axial force is transmitted over the length of the plate. The risk of a secondary loss of reduction is reduced.



Blood supply to the bone

Locking the screw into the plate does not generate plate-to-bone compression. Therefore, the periosteum will be protected and the blood supply to the bone preserved.



Combining Conventional and Locked Plating Techniques— The Locking Compression Plate (LCP) System


The combination of conventional compression plating and locked plating techniques enhances plate osteosynthesis. The result is a Combi hole that, depending on the indication, allows conventional compression plating, locked plating, or a combination of both.

Internal fixation using a combination of locking screws and standard screws

Note: If a combination of cortex and locking screws is used, a cortex screw should be inserted first to pull the plate to the bone.


If locking screws (1) have been used to fix a plate to a fragment, subsequent insertion of a conventional screw (2) in the same fragment without loosening and retightening the locking screw is not recommended.

Warning: If a locking screw is used first, care should be taken to ensure that the plate is held securely to the bone to avoid spinning of the plate about the bone.



Dynamic compression

In this example, once the metaphyseal fragment has been fixed with locking screws, the fracture can be dynamically compressed using a conventional screw in the DCU portion of the Combi hole.



Locked and conventional plating techniques

- Lag screws can be used to anatomically reduce the fracture and promote absolute stability.
- The behavior of a locking screw is not the same as that of a lag screw. With the locked plating technique, the implant locks the bone segments in their relative positions regardless of how they are reduced. Therefore, anatomical reduction must be achieved prior to implanting any locking screws.
- A plate used as a locked plate does not produce any additional compression between the plate and the bone.
- The unicortical insertion of locking screws creates a construct that is at least as strong as a construct made with bicortical insertion of conventional screws.

Screw Selection

The 4.0 mm and 5.0 mm locking screws are suitable for both diaphyseal and metaphyseal fractures. The 5.0 mm locking screw was designed as the principle screw for use with LCP Plates. It provides greater bending and shear strength than 4.0 mm locking screws (Chart 1). The 4.0 mm locking screw, with a 3.4 mm core diameter versus the 4.4 mm core diameter of the 5.0 mm locking screw, was developed to provide the option of placing a smaller diameter screw in small-statured patients or in cases where it is desirable to leave a smaller hole on explantation.

Locking screw fixation provides the greatest advantage over conventional screw fixation in poor quality bone. Even though 5.0 mm locking screws are significantly stronger in bending and shear than 4.0 mm locking screws, the behavior of both locking screw constructs provides relatively similar results in mechanical tests using 15 lb/ft³ foam, which simulates osteopenic bone, under axial load (Chart 2). Both bicortical locking screw constructs outperform a construct with conventional 4.5 mm cortex screws. When all constructs are tested in 40 lb/ft³ foam simulating good quality cortical bone, both locking and conventional constructs yield similar results when loaded axially (Chart 3).

Chart 1*




Chart 2*





Chart 3*



*DePuy Synthes. Data on File. Technical report: MT01-258: Static Testing of Locked Plates and Screws in Foam Material Locking Compression Plate Project Ref: P99-016. September 2001. Data represents test results from stainless steel implants only.

Preparation and Reduction

1. Plate selection

The plates are available in various lengths and configurations similar to the DePuy Synthes basic plate set. If necessary, use a bending template to determine plate length and configuration.

2. Contouring

Use the current bending instruments to contour the locking compression plate to the anatomy.

Note: The plate holes have been designed to accept some degree of deformation. When bending the plate, place the bending irons on two consecutive holes. This ensures that the threaded holes will not be distorted. Significant distortion of the locking holes will reduce locking effectiveness.

Precautions:

- **Reverse bending or use of the incorrect instrumentation for bending may weaken the plate and lead to premature plate failure (e.g. breakage). Do not bend the plate beyond what is required to match the anatomy.**
- **Do not bend the plate at the level of the holes.**

Note: Please refer also to the *AO Principles of Fracture Management*.¹

Precautions:

- **Instruments and screws may have sharp edges or moving joints that may pinch or tear user's glove or skin.**
- **Handle devices with care and dispose of worn bone cutting instruments in an approved sharps container.**

1. Rüedi TP, et al, ed., *AO Principles of Fracture Management*, New York: Thieme. 2000.

3. Reduction and temporary plate placement

Instruments

311.449 Push-Pull Reduction Device

324.075 Threaded Plate Holder

The plate may be temporarily held in place with standard plate holding forceps or the push-pull reduction device.

Note: The middle of the plate should be positioned over the fracture site if compression of the fracture fragments is desired.

Warning: For LCP Metaphyseal Plates, do not position the thinner portion of the plate over the fracture site.

The push-pull reduction device is designed to temporarily hold the plate to the bone through a plate hole. The device is self-drilling and connects with the DePuy Synthes quick coupling for power insertion. Insert into the near cortex only. After power insertion, turn the collet clockwise until it pulls the plate securely to the bone.

Note: Care should be taken to avoid inserting this device in a Combi hole that will be needed immediately for plate fixation. Also, if the chosen Combi hole is needed for placement of a locking screw, it is desirable to place the push-pull reduction device in the conventional portion of the Combi hole so that it does not interfere with the correct placement of the locking screw.

Alternatively, the threaded plate holder can be used as an aid to position the plate on the bone. The plate holder may also function as an insertion handle for use with minimally invasive plating techniques.



Screw Insertion

4. Screw insertion

Instrument

323.46	4.5 mm Universal Drill Guide
--------	------------------------------

Determine whether conventional cortex screws, cancellous bone screws, or locking screws will be used for fixation. A combination of all may be used.

Note: If a combination of cortex, cancellous and locking screws is used, a conventional screw should be used first to pull the plate to the bone.

Warning: If a locking screw is used first, care should be taken to ensure that the plate is held securely to the bone to avoid spinning of the plate about the bone as the locking screw is tightened to the plate.

Insertion of a cortex or cancellous bone screw

Use the 4.5 mm universal drill guide for an eccentric (compression) or neutral (buttress) insertion of cortex screws.

Note: The 4.5 mm LC-DCP Drill Guide (323.45) and the 4.5 mm DCP Drill Guide (322.44) are NOT suitable for use with LCP Plates.

Neutral insertion of a conventional screw

When pressing the universal drill guide into the DCU portion of the Combi hole, it will center itself and allow neutral predrilling.

Dynamic compression, eccentric insertion of a cortex screw

To drill a hole for dynamic compression, place the universal drill guide eccentrically at the edge of the DCU portion of the Combi hole, without applying pressure. Tightening of the cortex screws will result in dynamic compression corresponding to that of LC-DCP Plates.



Instruments

310.31	3.2 mm Drill Bit
310.431	4.3 mm Drill Bit
312.445	3.2 mm Threaded Drill Guide (for 4.0 mm screws)
312.449	4.3 mm Threaded Drill Guide (for 5.0 mm screws)
314.118	StarDrive Screwdriver, T25
314.119	StarDrive Screwdriver Shaft, T25, quick coupling
319.10	Depth Gauge
511.771*	Torque Limiting Attachment, 4 Nm or
511.774	Torque Limiting Attachment, 4 Nm, for AO Reaming Coupler

Insertion of 4.0 mm and 5.0 mm locking screws

Note: The locking screw is not a lag screw. Use nonlocking screws when requiring a precise anatomical reduction (e.g., joint surfaces) or interfragmentary compression. Before inserting the first locking screw, perform anatomical reduction and fix the fracture with lag screws, if necessary. After the insertion of locking screws, an anatomical reduction will no longer be possible without loosening the locking screw.

Screw the appropriate threaded drill guide for 4.0 mm screws and for 5.0 mm screws into an LCP Plate hole until fully seated (a). The use of a threaded drill guide is critical to ensure proper mating of the locking screw in the threaded portion of the Combi hole. The drill guide also has internal threads so guides can be assembled in series to increase length for percutaneous use.

Precaution: The use of the threaded drill guide is mandatory in order to ensure that the locking screw is drilled in the proper perpendicular angle and correctly locked in the plate.

Warning: Do not try to bend the plate using the threaded drill guide because damage may occur to the threads.

*Also available.



a

4. Screw insertion continued

Notes:

- Since the direction of a locking screw is determined by plate design, final screw position may be verified with a guide wire before insertion. This is especially important when the plate has been contoured or applied in metaphyseal regions around joint surfaces.
- The 5.0 mm cannulated locking screws and 5.0 mm cannulated conical screws for the locking periarticular plating system are compatible with the large fragment LCP Plates.

Use the appropriate drill bit (3.2 mm for 4.0 mm screws and 4.3 mm for 5.0 mm screws) to drill to the desired depth (b).

Remove the drill guide.

Precautions:

- Always irrigate during drilling to avoid thermal damage to the bone.
- For long screws and thick cortical bone, ensure sufficient cooling during insertion.
- Do not insert the drill tip into the opposite cortical bone to prevent damage of the opposite side structure and to avoid removal problems.

Use the depth gauge to determine screw length (c).

Insert the locking screw under power using the torque limiting attachment and StarDrive Screwdriver shaft (d).

The torque limiting attachment controls the tightening torque to 4 Nm.

- Ensures that enough torque is used to minimize the risk of the locking screw backing out of the plate
- Avoids locking the screw to the plate at full speed and minimizes the risk of cold-welding the screw to the plate
- DO NOT fully insert the locking screws by power without using the torque limiting attachment

Notes:

- The screw is securely locked to the plate when a click is heard.
- Do not lock the screws at full speed to reduce the risk of stripping the head. This can make it difficult to remove the implant.



Warning: Locking screws may be partially inserted using power equipment alone. However, never use power equipment to seat the locking screws into the plate without a torque limiting attachment (TLA).

Alternative method of locking screw insertion

Use the StarDrive Screwdriver to manually insert the appropriate locking screw. Carefully tighten the locking screw, as excessive force is not necessary to produce effective screw-to-plate locking.

Screw Placement Verification

5. Screw placement verification

Instruments

292.656	2.0 mm Non-Threaded Guide Wire
323.021	Direct Measuring Device
323.046	2.0 mm Wire Sleeve

Since the direction of a locking screw is determined by plate design, final screw position may be verified with a guide wire prior to insertion. This becomes especially important when the plate has been contoured or applied in metaphyseal regions around joint surfaces.

With the threaded drill guide in place, thread the 2.0 mm wire sleeve into the threaded drill guide until fully seated (a).

Insert a 2.0 mm nonthreaded guide wire through the wire sleeve to the desired depth (b).

- Verify guide wire placement under image intensification to determine if final screw placement is acceptable (c).

Notes:

- The guide wire position represents the final position of the locking screw. Confirm that the guide wire does not enter the joint.**
- If the angle of the locking screw is not optimal, it can be corrected. Bend the plate as needed, or move it in a proximal or distal direction. This technique is also suitable to preliminarily fix the plate to the bone.**

Measurement may be taken by sliding the tapered end of the direct measuring device over the guide wire down to the wire sleeve (d).

Remove the direct measuring device, guide wire and wire sleeve, leaving the threaded drill guide intact.

Use the appropriate size drill bit to drill the near cortex. Remove the threaded drill guide. Insert the appropriate length locking screw.



a



b



c



d

Postoperative Treatment and Implant Removal

Postoperative treatment

Postoperative treatment with locking compression plates does not differ from conventional internal fixation procedures.

Implant removal

To remove locking screws, unlock all screws from the plate; then remove the screws completely from the bone. This prevents simultaneous rotation of the plate when removing the last locking screw.

Screws

4.0 mm Locking Screws, self-tapping


Create a locked, fixed-angle screw/plate construct

- Available in 14 mm–62 mm lengths (2 mm increments)
- Threaded conical head
- Fully threaded shaft
- Self-tapping tip
- Available in stainless steel



4.5 mm Cortex Screws


- May be used in the DCU portion of the Combi holes in the plate shaft
- Compress the plate to the bone or create axial compression
- Self-tapping tip
- Available in stainless steel or titanium



5.0 mm Locking Screws, self-tapping

Create a locked, fixed-angle screw/plate construct

- Available in 14 mm–50 mm lengths (2 mm increments)
- Available in 55 mm–90 mm lengths (5 mm increments)
- Threaded conical head
- Fully threaded shaft
- Self-tapping tip
- Available in stainless steel or titanium alloy**



6.5 mm Cancellous Bone Screws, partially threaded

- Compress the plate to the bone and provide interfragmentary compression
- May be used in the DCU portion of the first and last Combi holes of the 4.5 mm narrow and broad LCP Plates
- Available in 16 mm, and 32 mm thread lengths
- Available in stainless steel or titanium



6.5 mm Cancellous Bone Screws, fully threaded

- Compress the plate to the bone
- May be used in the DCU portion of the first and last Combi holes of the 4.5 mm narrow and broad LCP Plates
- Available in stainless steel or titanium



**Ti-6Al-7Nb.

Plates

4.5 mm Narrow LCP Plates

- Available with 2–22 holes
- Available in stainless steel or titanium



224.591

4.5 mm Broad LCP Plates

- Available with 6–22 holes
- Available in stainless steel or titanium



226.591

4.5 mm LCP Proximal Tibia Plates*

- Available in stainless steel with 4, 6, 8, 10, 12, 14, 16, 18, and 20 shaft holes†
- Available in titanium with 4, 6, 8, 10, 12, and 14 shaft holes
- Available in left and right configurations



240.039

4.5 mm LCP T-Plates

- Available with 4, 6, and 8 shaft holes
- Available in stainless steel or titanium



240.161

*Also available.

†Plates with 16, 18, and 20 holes available sterile only.

Selected Instruments

310.31 3.2 mm Drill Bit



310.431 4.3 mm Drill Bit



312.445 3.2 mm Threaded Drill Guide



312.449 4.3 mm Threaded Drill Guide



314.118 StarDrive Screwdriver, T25



314.119 StarDrive Screwdriver Shaft, T25



Selected Instruments

323.021 Direct Measuring Device



323.046 2.0 mm Wire Sleeve



324.075 Threaded Plate Holder



397.706 Handle, for AO Reaming Coupler Connection



511.774 Torque Limiting Attachment (TLA), 4 Nm,
for AO Reaming Coupler




Large Fragment LCP Instrument and Implant Sets with 4.0 mm and 5.0 mm Locking Screws

Stainless Steel (115.400) and Titanium (146.400)


Graphic Case

- 690.360 Large Fragment LCP Implant Set Graphic Case (for stainless steel implants)
- 690.362 Large Fragment LCP Screw Set Graphic Case
- 690.363 Large Fragment LCP Instrument Set Graphic Case
- 690.420 Large Fragment LCP Titanium Implant Set Graphic Case



Instruments

- 292.656 2.0 mm Non-Threaded Guide Wire, 230 mm, spade point, 10 each
- 310.19 2.0 mm Drill Bit, 100 mm, quick coupling, 2 each
- 310.31 3.2 mm Drill Bit, 145 mm, quick coupling, 2 each
- 310.431* 4.3 mm Drill Bit, 180 mm, quick coupling, for 5.0 mm Locking Screws, 2 each
- 310.44 4.5 mm Drill Bit, 145 mm, quick coupling, 2 each
- 310.99 Countersink, for 4.5 mm and 6.5 mm screws
- 311.44 T-Handle, with quick coupling
- 311.449* Push-Pull Reduction Device, for use with 4.5 mm LCP plates, 2 each
- 311.46 Tap for 4.5 mm screws, 2 each
- 311.66 Tap for 6.5 mm Cancellous Bone Screws, 2 each
- 312.445* 3.2 mm Threaded Drill Guide, 2 each
- 312.449* 4.3 mm Threaded Drill Guide, 4 each
- 312.46 4.5 mm/3.2 mm Double Drill Sleeve
- 312.48 4.5 mm/3.2 mm Insert Drill Sleeve
- 312.67 6.5 mm/3.2 mm Double Drill Sleeve
- 314.11 Holding Sleeve
- 314.118* StarDrive Screwdriver, T25
- 314.119* StarDrive Screwdriver Shaft, T25, quick coupling
- 314.15 Large Hexagonal Screwdriver Shaft
- 314.27 Large Hexagonal Screwdriver
- 319.10 Depth Gauge, for 4.5 mm and 6.5 mm screws
- 319.39 Sharp Hook
- 319.97 Screw Forceps



*LCP System-specific instruments.

Note: For additional information, please refer to the package insert or www.e-ifu.com.

For detailed cleaning and sterilization instructions, please refer to www.depuysynthes.com/hcp/cleaning-sterilization or sterilization instructions, if provided in the instructions for use.

Large Fragment LCP Instrument and Implant Sets with 4.0 mm and 5.0 mm Locking Screws
Stainless Steel (115.400) and Titanium (146.400)

Instruments continued				5.0 mm Locking Screws, self-tapping, with T25 StarDrive Recess			
				Stainless			
				Steel	Titanium	Length (mm)	Quantity
321.12	Articulated Tensioning Device			212.201	412.201	14	5
321.15	Socket Wrench with universal joint, 11.0 mm width across flats			212.202	412.202	16	5
323.021*	Direct Measuring Device			212.203	412.203	18	5
323.046*	2.0 mm Wire Sleeve, 2 each			212.204	412.204	20	5
323.46	4.5 mm Universal Drill Guide			212.205	412.205	22	5
324.075*	Threaded Plate Holder			212.206	412.206	24	5
	Bending Templates			212.207	412.207	26	5
329.92	12 holes			212.208	412.208	28	5
329.97	7 holes			212.209	412.209	30	5
329.99	9 holes			212.210	412.210	32	5
397.706*	Handle, for AO Reaming Coupler connection			212.211	412.211	34	5
511.774*	Torque Limiting Attachment, 4 Nm, for AO Reaming Coupler			212.212	412.212	36	5
				212.213	412.213	38	5
				212.214	412.214	40	5
				212.215	412.215	42	5
				212.216	412.216	44	3
				212.217	412.217	46	3
				212.218	412.218	48	3
				212.219	412.219	50	3
				212.220	412.220	55	3
				212.221	412.221	60	3
				212.222	412.222	65	3
				212.223	412.223	70	3
				212.224	412.224	75	3
				212.225	412.225	80	3
				212.226	412.226	85	3
				212.227	412.227	90	3
Implants							
4.0 mm Locking Screws, self-tapping, with T25 StarDrive Recess, 3 each							
Stainless Steel		Length (mm)		Stainless Steel		Length (mm)	
02.204.014	14	02.204.040	40				
02.204.016	16	02.204.042	42				
02.204.018	18	02.204.044	44				
02.204.020	20	02.204.046	46				
02.204.022	22	02.204.048	48				
02.204.024	24	02.204.050	50				
02.204.026	26	02.204.052	52				
02.204.028	28	02.204.054	54				
02.204.030	30	02.204.056	56				
02.204.032	32	02.204.058	58				
02.204.034	34	02.204.060	60				
02.204.036	36	02.204.062	62				
02.204.038	38						

*LCP System-specific instruments.

Large Fragment LCP Instrument and Implant Sets with 4.0 mm and 5.0 mm Locking Screws
Stainless Steel (115.400) and Titanium (146.400)

Implants continued

4.5 mm Cortex Screws, self-tapping

Stainless

Steel	Titanium	Length (mm)	Quantity
214.814	414.814	14	4
214.816	414.816	16	4
214.818	414.818	18	4
214.820	414.820	20	6
214.822	414.822	22	6
214.824	414.824	24	6
214.826	414.826	26	12
214.828	414.828	28	12
214.830	414.830	30	12
214.832	414.832	32	12
214.834	414.834	34	12
214.836	414.836	36	12
214.838	414.838	38	12
214.840	414.840	40	12
214.842	414.842	42	12
214.844	414.844	44	4
214.846	414.846	46	4
214.848	414.848	48	4
214.850	414.850	50	4
214.852	414.852	52	4
214.854	414.854	54	4
214.856	414.856	56	4
214.858	414.858	58	4
214.860	414.860	60	4
214.862	414.862	62	4
214.864	414.864	64	4
214.866	414.866	66	4
214.868	414.868	68	4
214.870	414.870	70	4

6.5 mm Cancellous Bone Screws, 16 mm thread length, 3 each
Stainless

Steel	Titanium	Length (mm)
216.030	416.030	30
216.035	416.035	35
216.040	416.040	40
216.045	416.045	45
216.050	416.050	50
216.055	416.055	55
216.060	416.060	60
216.065	416.065	65
216.070	416.070	70
216.075	416.075	75
216.080	416.080	80
216.085	416.085	85
216.090	416.090	90
216.095	416.095	95
216.100	416.100	100
216.105	416.105	105
216.110	416.110	110

6.5 mm Cancellous Bone Screws, 32 mm thread length,
3 each

Stainless

Steel	Titanium	Length
217.045	417.045	45
217.050	417.050	50
217.055	417.055	55
217.060	417.060	60
217.065	417.065	65
217.070	417.070	70
217.075	417.075	75
217.080	417.080	80
217.085	417.085	85
217.090	417.090	90
217.095	417.095	95
217.100	417.100	100
217.105	417.105	105
217.110	417.110	110

Large Fragment LCP Instrument and Implant Sets with 4.0 mm and 5.0 mm Locking Screws
Stainless Steel (115.400) and Titanium (146.400)

Implants continued

6.5 mm Cancellous Bone Screws, fully threaded, 3 each

Stainless

Steel	Titanium	Length (mm)
218.025	418.025	25
218.030	418.030	30
218.035	418.035	35
218.040	418.040	40
218.045	418.045	45
218.050	418.050	50
218.055	418.055	55
218.060	418.060	60

Washer

Stainless

Steel	Titanium	Length (mm)	Quantity
219.99	419.99	13	6

4.5 mm Narrow LCP Plates

Stainless

Steel	Titanium	Holes	Length (mm)	Quantity
224.541	424.541	4	80	1
224.551	424.551	5	98	1
224.561	424.561	6	116	2
224.571	424.571	7	134	1
224.581	424.581	8	152	2
224.591	424.591	9	170	1
224.601	424.601	10	188	2
224.611	424.611	11	206	1
224.621	424.621	12	224	2
224.641	424.641	14	260	1
224.661	424.661	16	296	1

4.5 mm Broad LCP Plates

Stainless

Steel	Titanium	Holes	Length (mm)
226.561	426.561	6	116
226.571	426.571	7	134
226.581	426.581	8	152
226.591	426.591	9	170
226.601	426.601	10	188
226.611	426.611	11	206
226.621	426.621	12	224
226.641	426.641	14	260
226.661	426.661	16	296

4.5 mm LCP Medial Proximal Tibia Plates

Stainless

Steel	Titanium	Holes	Length (mm)
239.984	—	4, right	106
239.985	—	4, left	106
239.986	—	6, right	142
239.987	—	6, left	142

4.5 mm LCP T-Plates

Stainless

Steel	Titanium	Holes	Length (mm)
240.141	440.141	4	83
240.161	440.161	6	115
240.181	440.181	8	147

Also Available

Sets

01.225.602	LCP Periprosthetic System
01.225.604	Titanium LCP Periprosthetic System
01.106.902	LCP Metaphyseal Plate Set
01.106.904	Titanium LCP Metaphyseal Plate Set
01.120.432	4.5 mm LCP Medial Proximal Tibia Plate Implant Set
01.240.201	Periarticular LCP Plating System, with 5.0 mm Locking Screws
105.221	4.5 mm LCP Condylar Plate Implant Set
105.222	4.5 mm LCP Proximal Tibia Plate Implant Set
105.272	4.5 mm LCP Proximal Femur Plate Set
105.273	4.5 mm LCP Proximal Femur Hook Plate Set
115.401	Large Fragment LCP Instrument Set
115.402	Large Fragment LCP Plate Set
146.402	Large Fragment Titanium LCP Plate Set
115.403	Large Fragment LCP Screw Set
146.403	Large Fragment Titanium LCP Screw Set

Instruments

292.652	2.0 mm Threaded Guide Wire
397.705*	Handle, quick coupling, for ComPact Air Drive connection
511.771*	Torque Limiting Attachment, 4 Nm

Implants

02.001.300– 02.001.332	4.5 mm LCP Curved Condylar Plates, 10–22 holes, 242 mm–458 mm [◇]
222.656– 222.669	4.5 mm LCP Condylar Plates, 6–18 holes, 170 mm–386 mm [◇]

4.5 mm Narrow LCP Plates

Stainless			
Steel	Titanium	Holes	Length (mm)
224.521	424.521	2	44
224.531	424.531	3	62
224.631	424.631	13	242
224.651	424.651	15	278
224.681	424.681	18	332
224.701	424.701	20	368
224.721	424.721	22	404

224.753–
224.765 4.5 mm/3.5 mm LCP Metaphyseal Plates, 8–20 holes, 118 mm–334 mm

4.5 mm Curved Broad LCP Plates[◇]

Stainless		
Steel	Titanium	
226.622– 226.682	426.622– 426.682	12–18 holes, 229 mm–336 mm
226.702S– 226.762S	426.702S– 426.762S	20–26 holes [†] , 372 mm–479 mm

4.5 mm Broad LCP Plates

Stainless			
Steel	Titanium	Holes	Length (mm)
226.681	426.681	18	332
226.701	426.701	20	368
226.721	426.721	22	404

4.5 mm LCP Straight Reconstruction Plates

Stainless		
Steel	Titanium	
229.331– 229.461	429.331– 429.461	3–16 holes, 56 mm–303 mm

*LCP System-specific instruments.

[◇]These implants are available nonsterile and sterile-packed. Add "S" to catalog number for sterile product.

[†]Sterile only.

Implants continued

239.988– 4.5 mm LCP Medial Proximal Tibia Plates,
 239.997 8–16 holes, right and left,
 178 mm–322 mm[◇]

439.984– 4.5 mm Titanium LCP Medial Proximal Tibia
 439.997 Plates, 4–16 holes, right and left,
 106 mm–322 mm[◇]

4.5 mm LCP Proximal Tibia Plates[◇]

Stainless

Steel

Titanium

Steel	Titanium	
240.036–	440.036–	4–14 holes, right and left
240.047	440.047	82 mm–262 mm

240.0485–	—	16–20 holes, right and left [†]
240.0555		298 mm–370 mm

242.102– 4.5 mm LCP Proximal Femur Plates,
 242.116 2–16 holes, left, 139 mm–391 mm[◇]

242.802– 4.5 mm LCP Proximal Femur Plates,
 242.816 2–16 holes, right, 139 mm–391 mm[◇]

242.120– 4.5 mm LCP Proximal Femur Hook Plates,
 242.128 2–18 holes, 133 mm–421 mm

[◇] These implants are available nonsterile and sterile-packed. Add "S" to catalog number for sterile product.

[†] Sterile only.

Discontinued Product

Implants

4.0 mm Locking Screws, self-tapping, with T25 StarDrive Recess (Discontinued–December 2016)

Titanium

04.204.014
04.204.016
04.204.018
04.204.020
04.204.022
04.204.024
04.204.026
04.204.028
04.204.030
04.204.032
04.204.034
04.204.036
04.204.038
04.204.040
04.204.042
04.204.044
04.204.046
04.204.048
04.204.050
04.204.052
04.204.054
04.204.056
04.204.058
04.204.060
04.204.062

4.5 mm Shaft Screws (Discontinued–December 2016)

Stainless Steel

214.228
214.230
214.232
214.234
214.236
214.238
214.240
214.242
214.244
214.246
214.248
214.250

Titanium

414.228
414.230
414.232
414.234
414.236
414.238
414.240
414.242
414.244
414.246
414.248
414.250

Limited Warranty and Disclaimer: DePuy Synthes products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

Please also refer to the package insert(s) or other labeling associated with the devices identified in this surgical technique for additional information.

CAUTION: Federal Law restricts these devices to sale by or on the order of a physician.

Some devices listed in this technique guide may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.

Not all products may currently be available in all markets.



Manufactured or distributed by:

Synthes USA Products, LLC
1302 Wrights Lane East
West Chester, PA 19380

Synthes USA, LLC
1101 Synthes Avenue
Monument, CO 80132

Synthes GmbH
Luzernstrasse 21
4528 Zuchwil, Switzerland

To order (USA): 800-523-0322

To order (Canada): 844-243-4321

Note: For recognized manufacturer, refer to the product label.

www.depuysynthes.com