

TRUESPAN™ Meniscal Repair System Surgical Technique

Step 1: Remove from packaging and adjust white depth stop position as needed

The default setting for the depth stop is 18mm. It can be such that between 10mm and 20mm of the needle tip is exposed from the end of the tube. Depth stop position can be changed after the device has been inserted into the knee joint, and between implant deployments, if desired.



Step 2: Insert device into knee joint using a malleable graft retractor



Step 3: Penetrate tissue in desired location, press depth stop against tissue and pull trigger until click is heard, while maintaining forward pressure with the deployment gun



Step 4: Fully release trigger and remove device from tissue. Reposition and penetrate tissue with needle in desired location for 2nd implant, 6-9mm away from the 1st implant



Step 5: Pull trigger until click is heard to deploy 2nd implant. Release trigger and remove entire device from knee joint



Step 6: Pull suture with continuous force to tighten repair until desired tension is reached. Cut suture tail flush with tissue surface using a suture cutter. Repeat with additional devices as required

