

RIGIDLOOP™ CORTICAL FIXATION SYSTEM



15 mm Implant



XL Implant



Variable Depth Gauge



Disposable Kit

Designed with simplicity in mind:

- Offers enhanced strength while minimizing graft laxity¹
- Innovative variable depth gauge system eliminates tunnel and graft calculation complexities
- Convenience of a fully assembled procedural instrumentation kit
- Low profile reamers facilitate anatomic tunnel drilling through the medial portal while avoiding damage to the condyle


This publication is not intended for distribution in the USA.

SALES SHEET


RIGIDLOOP™
CORTICAL FIXATION SYSTEM

The RIGIDLOOP System offers enhanced strength and minimal elongation¹

Ultimate strength of loop construct (N)



Elongation of loop construct (mm)



RIGIDLOOP System ordering information

Implants


Cat No	Description
232009	RIGIDLOOP 15mm Fixed Loop Implant (Sterile)
232011	RIGIDLOOP 20mm Fixed Loop Implant (Sterile)
232014	RIGIDLOOP 25mm Fixed Loop Implant (Sterile)
232016	RIGIDLOOP 30mm Fixed Loop Implant (Sterile)
232017	RIGIDLOOP 35mm Fixed Loop Implant (Sterile)
232018	RIGIDLOOP 40mm Fixed Loop Implant (Sterile)
232019	RIGIDLOOP 45mm Fixed Loop Implant (Sterile)
232021	RIGIDLOOP 50mm Fixed Loop Implant (Sterile)
232022	RIGIDLOOP 55mm Fixed Loop Implant (Sterile)
232023	RIGIDLOOP 60mm Fixed Loop Implant (Sterile)
232034	RIGIDLOOP XL Button (Sterile)

Instruments

Cat No	Description
232037	RIGIDLOOP ACL Disposables Kit (Sterile)
232025	RIGIDLOOP Beath Pin (Non-Sterile)
232028	RIGIDLOOP 4.5mm Cannulated Reamer (Non-Sterile)
232030	RIGIDLOOP Standard Depth Gauge (Non-Sterile)
232031	RIGIDLOOP Variable Depth Gauge (Non-Sterile)
232032	RIGIDLOOP Instrument Caddy (Non-Sterile)
232033	RIGIDLOOP Silicon Mat for Instrument Caddy (Non-Sterile)
215451	Instrument Tray Lid (Non-Sterile)

Low Profile Reamers

Cat No	Description
219624	6.0mm Low Profile Reamer (Non-Sterile)
219625	6.5mm Low Profile Reamer (Non-Sterile)
219626	7.0mm Low Profile Reamer (Non-Sterile)
219627	7.5mm Low Profile Reamer (Non-Sterile)
219628	8.0mm Low Profile Reamer (Non-Sterile)
219629	8.5mm Low Profile Reamer (Non-Sterile)
219630	9.0mm Low Profile Reamer (Non-Sterile)
219631	9.5mm Low Profile Reamer (Non-Sterile)
219632	10.0mm Low Profile Reamer (Non-Sterile)
219633	10.5mm Low Profile Reamer (Non-Sterile)
219634	11.0mm Low Profile Reamer (Non-Sterile)
219635	11.5mm Low Profile Reamer (Non-Sterile)
219636	12.0mm Low Profile Reamer (Non-Sterile)
219637	13.0mm Low Profile Reamer (Non-Sterile)



Low Profile Reamer

References

1. Internal test data, ref. DHF-102035-MEMO14 &16.

Limited Warranty and Disclaimer: DePuy Mitek, Inc. products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

Not all products are currently available in all markets.

Third party trademarks used herein are trademarks of their respective owners.



Johnson & Johnson Medical Limited PO BOX 1988, Simpson Parkway, Livingston, West Lothian, EH54 0AB, United Kingdom. Incorporated and registered in Scotland under company number SC132162.

DePuy Mitek, Inc.
325 Paramount Drive
Raynham, MA 02767-0350
USA

Medos International SÀRL
Chemin-Blanc 38
2400 Le Locle
Switzerland

Distributed in the USA by:
DePuy Mitek, Inc.
325 Paramount Dr.
Raynham, MA 02767
USA
Tel: +1 (800) 382-4682

Authorized European Representative:
DePuy International Ltd
St Anthony's Road
Leeds LS11 8DT
England
Tel: +44 (0)113 270 0461
Fax: +44 (0)113 272 4101



depuysynthes.com

©Johnson & Johnson Medical Limited. 2016. All rights reserved.

CA#DSEM/MTK/0216/0305 Issued: 03/16