


ACRION™ SAW BLADES CROSS REFERENCE CHART

This chart shows the DePuy Synthes ACRION™ Saw Blades that have replaced the Joint and Stablecut® Saw Blades. The teeth geometry of Oscillating ACRION™ Saw Blades was developed by DePuy Synthes.



Oscillating Saw Blades with DePuy Synthes / Zimmer® Coupling

ACRION™ Saw Blades SKUs				Current SKUs			
Part number	Usable length [mm]	Width [mm]	Cutting thickness [mm]	Joint Saw Blade Part number		Stablecut Saw Blade Part number	
05.003.400S	60	13	0.89				
05.003.401S	90	13	0.89	05.003.100S			
05.003.402S	90	19	0.89	05.003.110S	05.003.120S		
05.003.403S	90	25	0.89	05.003.130S		9640-05-000	
05.003.411S	90	13	1.00	05.003.102S	05.003.103S		
05.003.413S	90	25	1.00	05.003.132S	05.003.133S		
05.003.421S	90	13	1.19	05.003.104S		5627-21-500	9640-01-820
05.003.422S	90	19	1.19	05.003.114S			
05.003.423S	90	25	1.19	05.003.134S		9640-01-826	5627-20-500
05.003.424S	90	13-8	1.19			5627-09-552	
05.003.425S	90	21-13	1.19	05.003.124S		5627-09-550	
05.003.431S	90	13	1.27	05.003.106S			
05.003.432S	90	19	1.27	05.003.116S			
05.003.433S	90	25	1.27	05.003.136S		5627-79-000	
05.003.434S	70	13-8	1.27			5627-90-000	
05.003.435S	90	21-13	1.27	05.003.126S			
05.003.442S	90	19	1.37	05.003.117S	05.003.127S		
05.003.443S	90	25	1.37	05.003.137S			
05.003.451S	90	13	1.47	05.003.109S			
05.003.452S	90	19	1.47	05.003.119S	05.003.129S		
05.003.453S	90	25	1.47	05.003.139S		5627-37-000	9640-27-000
05.003.454S	70	13-8	1.47			9505-01-595	

DePuy Synthes / Zimmer coupling:



For use with:

- Battery Power Line II Oscillator (530.710)
- Trauma Recon System Modular (05.001.201) and oscillating saw attachments (05.001.223* or 05.001.224)
- Trauma Recon System Recon Saw (05.001.240)
- Compact Air Drive II (511.701) and oscillating saw attachments (511.800, 511.80 US part number or 511.801)
- Zimmer Biomet: Zimmer Universal Power System (89-8509-450-60)

*05.001.223 not available in the USA.

Oscillating Saw Blades with Stryker® Coupling

ACRION™ Saw Blades SKUs				Current SKUs			
Part number	Usable length [mm]	Width [mm]	Cutting thickness [mm]	Stablecut Saw Blade Part number			
05.003.500S	60	13	0.89	5627-44-000	5627-44-601		
05.003.521S	90	13	1.19	5627-11-500	5627-11-601	9640-01-818	
05.003.523S	90	25	1.19	5627-10-500	5627-10-601	9640-01-823	9640-01-824
05.003.524S	90	13-8	1.19	5627-10-552	5627-10-652		
05.003.525S	90	21-13	1.19	5627-10-651	5627-10-550		
05.003.533S	90	25	1.27	5627-82-601	5627-82-000		
05.003.534S	70	13-8	1.27	5627-87-000	5627-87-601		
05.003.535S	90	21-13	1.27				
05.003.553S	90	25	1.47	5627-53-000	5627-53-601	9640-20-000	9640-23-000
05.003.554S	70	13-8	1.47	9505-01-591	9505-01-592	9505-01-924	

Stryker coupling:



For use with:

- System EHD (296-98), System 2000 (2108)
- System 4 (4108), System 5 (4208), System 6 (6208), System 7 (7208), System 8 (8208)

Oscillating Saw Blades with ConMed Linvatec® Hall® Coupling

ACRION™ Saw Blades SKUs				Current SKUs		
Part number	Usable length [mm]	Width [mm]	Cutting thickness [mm]	Stablecut Saw Blade Part number		
05.003.621S	90	13	1.19	5627-15-501	9640-01-816	9640-01-822
05.003.623S	90	25	1.19	5627-14-501	9640-01-827	
05.003.624S	90	13-8	1.19	5627-14-552		
05.003.625S	90	21-13	1.19	5627-14-550	5627-15-550	
05.003.633S	90	25	1.27	5627-84-501		
05.003.634S	90	13-8	1.27			
05.003.635S	90	21-13	1.27			
05.003.653S	90	25	1.47	5627-55-501	9640-21-100	9640-28-000
05.003.654S	70	13-8	1.47	9505-01-593	9505-01-594	

ConMed Linvatec® Hall® coupling:



For use with:

- MPower® (PRO6300), MPower® 2 (PRO6350M)
- PowerPro Max™ (PRO5300M), PowerPro® & PowerPro® Electric II (PRO6125)
- PowerPro Pneumatic (PRO6175), Hall Series 3 (5044-002), Hall Series 4 (5067-002)
- Versipower Plus (5071-002), Hall 50™ (PRO7300B)

Oscillating Saw Blades with Aesculap® Coupling

ACRION™ Saw Blades SKUs			
Part number	Usable length [mm]	Width [mm]	Cutting thickness [mm]
05.003.703S	90	25	0.89
05.003.721S	90	13	1.19
05.003.723S	90	25	1.19
05.003.733S	90	25	1.27
05.003.734S	90	13-8	1.27
05.003.735S	90	21-13	1.27
05.003.753S	90	25	1.47

Aesculap coupling:



For use with:

- Acculan 3Ti (GA673, GA668)
- Acculan 4 (GA331, GA336)

Aesculap saw blades are not available in USA.

Reciprocating Saw Blades

ACRION™ Saw Blades SKUs				Current SKUs	
Part number	Usable length [mm]	Width [mm]	Cutting thickness [mm]	Saw Blade Coupling	Stablecut Saw Blade Part number
511.917S	76	12.5	1.19	DePuy Synthes	
511.918S	76	12.5	1.19	Stryker	
511.919S	76	12.5	1.19	Linvatec Hall	
511.920S	76	12.5	1.47	DePuy Synthes	
511.921S	76	12.5	1.47	Stryker	2299-10-002
511.922S	76	12.5	1.47	Linvatec Hall	2299-10-001

DePuy Synthes coupling for use with:


- Battery Power Line II Reciprocator (530.715)
- Trauma Recon System Modular (05.001.201) and reciprocating saw attachment (05.001.225)
- Compact Air Drive II (511.701) and reciprocating saw attachment (511.902)

Stryker® coupling for use with:

- System EHD (296-96), System 2000 (2106)
- System 4 (4106), System 5 (4206), System 6 (6206), System 7 (7206), System 8 (8206)

ConMed Linvatec® Hall® coupling for use with:

- Hall 50™ (PRO7400B), MPower® (PRO6400)
- MPower® 2 (PRO6400M), PowerPro Pneumatic (PRO6185)
- PowerPro® Electric II (PRO6135), PowerPro Max™ (PRO5400M)
- Hall Series 4 (5067-003)



Limited Warranty and Disclaimer: DePuy Synthes products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

Please see the package insert(s), surgical technique guide(s), or other product labeling associated with the devices contained in this brochure for information on indications, warnings, precautions, and instructions for use.

Some devices listed in this brochure may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultants for items approved for sale in Canada.

Not all products are currently available in all markets. Data available on request.

All saw blades are supplied exclusively in sterile packaging as single use saw blades.

Detailed information on the ACRION Saw Blades can be found in the Large Bone Saw Blades brochure. Prior to ordering saw blades please check the exact dimension. This reference chart is for guidance only and subject to change.

The legal manufacturer of the ACRION Saw Blades is Gebr. Brasseler GmbH & Co. KG Germany. The blades are exclusively distributed by Synthes GmbH and The Anspach Effort, Inc.

The third party trademarks used herein are trademarks of their respective owners.

Please contact your DePuy Synthes sales representative for more information.

For further information on DePuy Synthes products please visit: www.depuyssynthes.com