



ANSPACH™ Power Tools

G1 Dissection Tools


Reference Chart



ANSPACH DISSECTION TOOLS

Reference Chart

Combine the prefix and suffix from the chart below to determine the dissection tool Catalog Number.
Example: S-1B-G1



FLUTED BALL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
05B-G1	0.5 mm Fluted Ball, Extended Length	•											
1B-G1	1 mm Fluted Ball, Extended Length	•		•									
15B-G1	1 mm Fluted Ball, Standard Length	•					•						
15B-G1	1.5 mm Fluted Ball, Extended Length	•											
155B-G1	1.5 mm Fluted Ball, Standard Length	•											
2B-G1	2 mm Fluted Ball, Extended Length	•		•									
25B-G1	2 mm Fluted Ball, Standard Length	•					•						
25BT-G1	2 mm Fluted Ball, Standard Length, Tapered	•											
25B-G1	2.5 mm Fluted Ball, Extended Length	•											
255B-G1	2.5 mm Fluted Ball, Standard Length	•											
35B-2-G1	3 mm Fluted Ball, Standard Length	•		•									
3B-G1	3 mm Fluted Ball, Extended Length	•		•									
35B-G1	3 mm Fluted Ball, Standard Length	•		•			•	•	•	•			
35BT-G1	3 mm Fluted Ball, Standard Length, Tapered	•											
35B-G1	3.5 mm Fluted Ball, Extended Length	•											
355B-G1	3.5 mm Fluted Ball, Standard Length	•											
4B-27-G1	4 mm Fluted Ball, Extended Length	•		•									
4BT-G1	4 mm Fluted Ball, Standard Length, Tapered	•											
4B-G1	4 mm Fluted Ball, Standard Length	•		•			•	•	•	•			
5B-12F-G1	5 mm Fluted Ball, 12 Flutes	•		•									
5B-2-G1	5 mm Fluted Ball, Standard Length	•		•									
5B-27-G1	5 mm Fluted Ball, Extended Length	•											
5B-A-G1	5 mm Fluted Ball, Less Aggressive	•											
5B-G1	5 mm Fluted Ball, Standard Length	•	•	•			•	•	•	•			
6B-G1	6 mm Fluted Ball, Standard Length	•	•	•	•		•	•	•	•			
7B-G1	7 mm Fluted Ball, Standard Length	•	•	•	•	•		•	•	•			
8B-G1	8 mm Fluted Ball, Standard Length	•											
9B-G1	9 mm Fluted Ball, Standard Length	•			•								

ANSPACH HELIX™ FLUTED BALL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
25B-HX-G1	ANSPACH HELIX™ 2 mm Fluted Ball, Standard Length	•		•									
2B-HX-G1	ANSPACH HELIX™ 2 mm Fluted Ball, Extended Length	•		•									
255B-HX-G1	ANSPACH HELIX™ 2.5 mm Fluted Ball, Standard Length	•											
25B-HX-G1	ANSPACH HELIX™ 2.5 mm Fluted Ball, Extended Length	•											
35B-HX-G1	ANSPACH HELIX™ 3 mm Fluted Ball, Standard Length	•		•									
3B-HX-G1	ANSPACH HELIX™ 3 mm Fluted Ball, Extended Length	•		•									
355B-HX-G1	ANSPACH HELIX™ 3.5 mm Fluted Ball, Standard Length	•											
35B-HX-G1	ANSPACH HELIX™ 3.5 mm Fluted Ball, Extended Length	•											
4B-HX-G1	ANSPACH HELIX™ 4 mm Fluted Ball, Standard Length	•		•									
5B-HX-G1	ANSPACH HELIX™ 5 mm Fluted Ball, Standard Length	•		•									
6B-HX-G1	ANSPACH HELIX™ 6 mm Fluted Ball, Standard Length	•		•									
7B-HX-G1	ANSPACH HELIX™ 7 mm Fluted Ball, Standard Length	•		•									

CARBIDE FLUTED BALL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
1B-C-G1	1 mm Carbide Fluted Ball, Extended Length	•											
15B-C-G1	1 mm Carbide Fluted Ball, Standard Length	•											
2B-C-G1	2 mm Carbide Fluted Ball, Extended Length	•											
25B-C-G1	2 mm Carbide Fluted Ball, Standard Length	•											
25B-C-BF-G1	2 mm Carbide Fluted Ball, Standard, 8 Flutes	•											
3B-C-G1	3 mm Carbide Fluted Ball, Extended Length	•											
35B-C-G1	3 mm Carbide Fluted Ball, Standard Length	•											
4B-C-10F-G1	4 mm Carbide Fluted Ball, 10 Flutes	•											
4B-C-G1	4 mm Carbide Fluted Ball	•											
4B-C-G1	4 mm Carbide Fluted Ball	•											
5B-C-10F-G1	5 mm Carbide Fluted Ball, 10 Flutes	•											
5B-C-G1	5 mm Carbide Fluted Ball	•											
6B-C-G1	6 mm Carbide Fluted Ball	•											
7B-C-14F-G1	7 mm Carbide Fluted Ball, 14 Flutes	•											
7B-C-G1	7 mm Carbide Fluted Ball	•											

DIAMOND BALL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
05D-G1	0.5 mm Diamond Ball, Extended Length	•											
055D-G1	0.5 mm Diamond Ball, Standard Length	•											
065D-G1	0.6 mm Diamond Ball, Standard Length	•											
085D-G1	0.8 mm Diamond Ball, Standard Length	•											
1D-G1	1 mm Diamond Ball, Extended Length	•		•									
15D-G1	1 mm Diamond Ball, Standard Length	•					•						
12D-G1	1.2 mm Diamond Ball, Extended Length	•											
15D-G1	1.5 mm Diamond Ball, Extended Length	•		•									
155D-G1	1.5 mm Diamond Ball, Standard Length	•											
25D-120-G1	2 mm Diamond Ball, 120/140 Grit, Standard Length	•											
2D-G1	2 mm Diamond Ball, Extended Length	•		•									
25D-G1	2 mm Diamond Ball, Standard Length	•						•					
25D-G1	2.5 mm Diamond Ball, Extended Length	•											
255D-G1	2.5 mm Diamond Ball, Standard Length	•											
35D-100-G1	3 mm Diamond Ball, 100/120 Grit, Standard Length	•											
3D-G1	3 mm Diamond Ball, Extended Length	•											
35D-G1	3 mm Diamond Ball, Standard Length	•						•					
35D-G1	3.5 mm Diamond Ball, Extended Length	•											
355D-G1	3.5 mm Diamond Ball, Standard Length	•											
4D-100-G1	4 mm Diamond Ball, 100/120 Grit	•											
4D-L-G1	4 mm Diamond Ball	•											
4D-G1	4 mm Diamond Ball	•		•	•	•							
5D-100-G1	5 mm Diamond Ball, 100/120 Grit	•											
5D-2-G1	5 mm Diamond Ball	•											
5D-G1	5 mm Diamond Ball	•		•	•	•							
6D-100-G1	6 mm Diamond Ball, 100/120 Grit	•											
6D-G1	6 mm Diamond Ball	•											


COARSE DIAMOND BALL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
1DC-G1	1 mm Coarse Diamond Ball, Extended Length	•											
15DC-G1	1 mm Coarse Diamond Ball, Standard Length	•											
15DC-G1	1.5 mm Coarse Diamond Ball, Extended Length	•											
155DC-G1	1.5 mm Coarse Diamond Ball, Standard Length	•											
15DCT-G1	1.5 mm Coarse Diamond Ball, Tapered, Extended Length	•											
2DC-G1	2 mm Coarse Diamond Ball, Extended Length	•											
25DC-G1	2 mm Coarse Diamond Ball, Standard Length	•											
25DCT-G1	2 mm Coarse Diamond Ball, Tapered, Standard Length	•											
25DC-G1	2.5 mm Coarse Diamond Ball, Extended Length	•											
255DC-G1	2.5 mm Coarse Diamond Ball, Standard Length	•											
3DC-2-G1	3 mm Coarse Diamond Ball	•											
35DC-22-G1	3 mm Coarse Diamond Ball	•											
35DC-50-G1	3 mm Coarse Diamond Ball, 50/60 Grit, Standard Length	•											
3DC-G1	3 mm Coarse Diamond Ball, Extended Length	•											
35DC-G1	3 mm Coarse Diamond Ball, Standard Length	•						•					
35DCT-G1	3 mm Coarse Diamond Ball, Tapered, Standard Length	•											
35DC-G1	3.5 mm Coarse Diamond Ball, Extended Length	•											
355DC-G1	3.5 mm Coarse Diamond Ball, Standard Length	•											
4DC-22-G1	4 mm Coarse Diamond Ball	•											
4DC-G1	4 mm Coarse Diamond Ball	•							•				
5DC-22-G1	5 mm Coarse Diamond Ball	•											
5DC-40-G1	5 mm Coarse Diamond Ball, 40/60 Grit	•											
5DC-G1	5 mm Coarse Diamond Ball	•							•				
55DC-G1	5.5 mm Coarse Diamond Ball	•											
6DC-2-G1	6 mm Coarse Diamond Ball	•											
6DC-G1	6 mm Coarse Diamond Ball	•							•				
8DC-G1	8 mm Coarse Diamond Ball	•											

EXTRA COARSE DIAMOND BALL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
4DX-20-G1	4 mm Extra Coarse Diamond Ball, 20/30 Grit	•											
4DX-G1	4 mm Extra Coarse Diamond Ball	•											
5DX-G1	5 mm Extra Coarse Diamond Ball	•											
6DX-20-G1	6 mm Extra Coarse Diamond Ball, 20/30 Grit	•											
6DX-G1	6 mm Extra Coarse Diamond Ball	•											
7DX-G1	7 mm Extra Coarse Diamond Ball	•											
8DX-G1	8 mm Extra Coarse Diamond Ball	•											

Combine the prefix and suffix from the chart below to determine the dissection tool Catalog Number.
Example: S-1B-G1



ROUTER


Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
6R-G1	1.12 mm x 6.35 mm Fluted Router	•		•									
1R-G1	1.2 mm x 12.8 mm Fluted Router	•											
1R-G1	1.4 mm x 12.8 mm Fluted Router	•											
CRN-M-G1	1.3 mm x 12.8 mm Fluted Router	•									•		
CRN-M-G1	2 mm x 15.6 mm Fluted Router	•										•	
CRN-G1	1.4 mm x 12.8 mm Fluted Router	•										•	
CRN-G1	1.85 mm x 16 mm Fluted Router	•											•
CRN-G1	2.5 mm x 25.4 mm Fluted Router	•											•
14R-G1	1.4 mm x 12.8 mm Fluted Router	•			•								
CRN-S-G1	2.15 mm x 22 mm Fluted Spiral Router	•										•	
CRN-S-G1	2.47 mm x 30.3 mm Fluted Spiral Router	•											•
25R-G1	2.15 mm x 25.4 mm Fluted Router	•				•							
3-10-G1	1.93 mm Suture Pass Drill / 1.17 eye	•											

OVAL

Catalog Number	Description	S	SHD	M	MHD	L	LHD	MIA16	XLHD	ADG	P	A	L
40VDC-G1													

Minimal Access Attachments

Combine the prefix and suffix from the chart below to determine the dissection tool Catalog Number.
Example: MA10-1SB



Catalog Number Suffix **Catalog Number Prefix** →

ATTACHMENT TYPE

		MA10	MA15	MA19
FLUTED BALL				
1SB	1 mm Fluted Ball	●		
15SB	1.5 mm Fluted Ball	●		
2SB	2 mm Fluted Ball	●	●	●
25SB	2.5 mm Fluted Ball	●		
3SB	3 mm Fluted Ball	●	●	●
DIAMOND BALL				
05SD	0.5 mm Diamond Ball	●		
08SD	0.8 mm Diamond Ball	●		
15SD	1.5 mm Diamond Ball		●	●
2SD	2 mm Diamond Ball	●	●	●
3SD	3 mm Diamond Ball	●	●	●
4D	4 mm Diamond Ball	●	●	●
4DC	4 mm Coarse Diamond Ball		●	
MATCHSTICK				
8N-18	1.8 mm Fluted Matchstick	●	●	●
8N	2.3 mm Fluted Matchstick	●	●	●
8NS	3 mm Fluted Matchstick	●	●	●
8NSDC	3 mm Coarse Diamond Matchstick	●	●	●

OCM-G1 OTOLOGIC Curved Micro Attachment, G1

DIAMOND BALL CURVED BURRS

Burr Diameter	5.6 cm Sleeve Length	7.3 cm Sleeve Length
0.5 mm	OCM5-05SD	OCM7-05SD
0.6 mm	OCM5-06SD	OCM7-06SD
0.7 mm	OCM5-07SD	OCM7-07SD
0.8 mm	OCM5-08SD	OCM7-08SD
1.0 mm	OCM5-1SD	OCM7-1SD
1.5 mm	OCM5-15SD	OCM7-15SD
2.0 mm	OCM5-2SD	OCM7-2SD

FLUTED BALL CURVED BURRS

Burr Diameter	5.6 cm Sleeve Length	7.3 cm Sleeve Length
1.0 mm	OCM5-1SB	OCM7-1SB
1.5 mm	OCM5-15SB	OCM7-15SB
2.0 mm	OCM5-2SB	OCM7-2SB

For detailed information refer to the Instructions for Use of the corresponding Power Tools system.

Not all products are currently available in all markets.

This publication is not intended for distribution in the USA.

The third-party trademarks used herein are the trademarks of their respective owners.



The Anspach Effort, Inc.
4500 Riverside Drive
Palm Beach Gardens, FL 33410
USA

EC REP
Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland



www.jnjmedicaldevices.com