


Mitek Sports Medicine's All-Star Lineup

- 1 **SPEEDTRAP™** Graft Prep System
- 2 **TRUESPAN™** Meniscal Repair System
- 3 **TWISTR™** Retrograde Reamer & Cruciate+ Instrumentation
- 4 **INTRAFIX™ ADVANCE** Tibial Fastener System
- 5 **RIGIDLOOP™** Adjustable Cortical System



EFFICIENT. VERSATILE. REPRODUCIBLE.

Introducing the all-star team designed for streamlined and more efficient ACL & PCL reconstruction and meniscal repair procedures.

At DePuy Synthes Mitek Sports Medicine, we know that every second is valuable in your OR. That's why we are enhancing our commitment to streamlined and reproducible procedures. From graft preparation, meniscal repair, and retrograde reaming to femoral and tibial implants, Mitek Sports Medicine now offers an all-star lineup of solutions for your cruciate ligament reconstructions and meniscal repairs. Together with our comprehensive portfolio of precision instruments, we've got the solutions you need for a more efficient OR.



SPEEDTRAP™
Graft Prep System



TRUESPAN™
Meniscal Repair System



TWISTR™
Retrograde Reamer
& Cruciate+
Instrumentation



**INTRAFIX™
ADVANCE** Tibial
Fastener System



RIGIDLOOP™
Adjustable Cortical
System

Procedural assistance FROM START TO FINISH.

Explore our all-star team of products for cruciate ligament reconstruction and meniscal repair procedures and bring efficiency and reproducibility to your OR



SPEEDTRAP™ Graft Prep System:

The graft system designed to prepare grafts faster than traditional whipstitching techniques **by at least 77%**.¹



TRUESPAN™ Meniscal Repair System:

An all-inside meniscal repair device that features a **streamlined handle, single trigger, and auto-reloading mechanism** to minimize device 'fiddle-factor' and enable designed simple one-handed use.



TWISTR™ Retrograde Reamer & Cruciate+ Instrumentation:

An **adjustable, one-size-fits-all** retrograde reamer and instrumentation set for precision outside-in drilling during cruciate ligament reconstruction.



INTRAFIX™ ADVANCE Tibial Fastener System:

The sheath and screw implant system **designed for greater reliability, simplicity, and speed** while maintaining the biomechanical performance of the INTRAFIX™ ACL Tibial Fastener System.



RIGIDLOOP™ Adjustable Cortical System:

An adjustable cortical button that **combines the convenience and versatility of adjustable loop devices with rigid graft fixation comparable to Endobutton**.²

Strengthen your knee procedures with this all-star line up.

Call your Mitek Sports Medicine representative or visit www.depuysynthes.com today for more information.

References:

1. Reference Internal Data 103304691
2. Reference Internal Data 103356305