

# TRUESPAN™ MENISCAL REPAIR SYSTEM



**Meniscal Repair, *Simplified***

Getting your patients back to  
their passion is our passion.






**Simplified Deployment** – Device features an ergonomic design that is comfortable to hold and a single trigger, allowing use with only one hand! No needle assembly or reloading steps required.

**Robust** – Implants and suture are protected by a depth stop and rigid needle while accessing hard to reach locations.

**Versatile** – Available 0°, 12°, & 24° needles pre-loaded with non-absorbable PEEK or absorbable PLGA implants and partially absorbable 2-0 ORTHOCORD® Suture.

**Low Profile** – Small size needle & implants minimize tissue disruption and ‘kickback’ during deployment. Repairs with TRUESPAN™ Implants do not leave a knot or hard body on articular surfaces!



### Ordering Information

	Order Number	Product Information
PEEK IMPLANT SYSTEMS	228150	TRUESPAN™ Meniscal Repair System with PEEK Implants and 0° Needle
	228151	TRUESPAN™ Meniscal Repair System with PEEK Implants and 12° Needle
	228152	TRUESPAN™ Meniscal Repair System with PEEK Implants and 24° Needle
PLGA IMPLANT SYSTEMS	228160	TRUESPAN™ Meniscal Repair System with PLGA Implants and 0° Needle
	228161	TRUESPAN™ Meniscal Repair System with PLGA Implants and 12° Needle
	228162	TRUESPAN™ Meniscal Repair System with PLGA Implants and 24° Needle

	Order Number	Product Description
Related Products	232024	Malleable Graft Retractor, 5 Pack
	228302	Arthroscopic Pusher/Cutter
	253110	Meniscal Rasp Up 45°
	253111	Meniscal Rasp Up 90°
	253112	Meniscal Rasp Straight 90°, Calibrated
	253113	Meniscal Probe
	228001	Meniscal Depth Probe
	228144	#2/0 ORTHOCORD® with Double-arm Meniscal Needles
	228145	#2/0 ORTHOCORD® (three free strands) with Meniscal Needle
	214101	CHIA PERCPASSER® Suture Passer (5/box)

**Limited Warranty and Disclaimer:** DePuy Mitek, Inc. products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

**WARNING:** In the USA, this product has labeling limitations. See package insert for complete information.

**CAUTION:** USA Law restricts these devices to sale by or on the order of a physician.

To order in the US, call Johnson & Johnson Health Care Systems Inc. Customer Support Services at 800-255-2500.



**Medos International SÄRL**  
Chemin-Blanc 38  
2400 Le Locle  
Switzerland

[www.depuyssynthes.com](http://www.depuyssynthes.com)

*Distributed in the USA by:*  
**DePuy Mitek, Inc.**  
325 Paramount Drive  
Raynham, MA 02767 USA  
Tel: +1 (800) 382-4682